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ZALĄCZNIK A

Image 1 - Excessive Charge Time EOL (Observation)

OBSERVATIONS (5)

- ATP and shock therapies will not be delivered: charge circuit inactive. Inform a Medtronic rep.
- EOL: replace device immediately.
- Patient Alert: charge time was > 30 sec.
- Patient Alert: charge circuit timeout occurred.
- Patient Activity less than 2 hr/day for 2 weeks.

Image 2 – Excessive Charge Time EOL Alert (Programming Screen)

Patient Alert Setup

Sound tone for:

Lead Impedance Out of Range...	On-High
Low Battery Voltage ERI...	On-High
Excessive Charge Time EOL...	On-High
Number of Shocks Delivered in an Episode...	Off
All Therapies in a Zone Exhausted for an Episode...	Off
VF Detection OFF, 3+ VF or 3+ FVT Rx Off.	On-High
Alert Time...	09:10

Buttons: Demonstrate Tones..., Undo Pending, OK

Image 3 – Programming Steps to Change Auto-Cap Formation Interval to Fixed value of 6-month intervals

VF Therapies

	Rx1	Rx2	Rx3	Rx4	Rx5	Rx6
VF Therapy Status	On	On	On	On	On	On
Energy	30 J	30 J	30 J	30 J	30 J	30 J
Pathway	B>AX	B>AX	B>AX	B>AX	B>AX	AX>B
ATP...	During Charging					

Shared Settings... Auto Cap Formation...

Auto Cap Formation

Minimum Auto Cap Formation Interval: Auto

Note: Re-programming this parameter resets the Auto Cap Formation Interval timer.

Last Charge to Full Energy: 06-Jan-1994 04:42:11

Last Capacitor Formation: 06-Jan-1994 04:27:14

Buttons: Undo Pending, OK

Auto Cap Formation

Minimum Auto Cap Formation Interval: 6 month

Note: Re-programming this parameter resets the Auto Cap Formation Interval timer.

Auto Cap Formation Interval will be extended temporarily by 30 J charges.

Last Charge to Full Energy: 06-Jan-1994 04:42:11

Last Capacitor Formation: 06-Jan-1994 04:27:14

Buttons: Undo Pending, OK

Minimum Auto Cap Formation Interval

Auto [X] 2 4 6

1 3 5

- This value will not allow the device to make automatic cap formation adjustments.

Buttons: Undo Pending, Close

Images 4a and 4b – Programming Screens to Conduct In-clinic High-Voltage Test Charge

AAI+ VF VT Resume Suspend

72 bpm / 830 ms

ECG Lead II

EGM1: Atp/Aring

ECG2: Vtp/Vring

Tests - Charge/Dump

Last Charge to Full Energy: 06-Jan-1994 04:42:11

Charge Time: 4.6 sec

Energy: 0.1 - 30 J

Last Capacitor Formation*: 06-Jan-1994 04:27:14

Buttons: CHARGE, DUMP, Capactors, Capactors, Print...

AAI+ VF VT Resume Suspend Manual operation charging

72 bpm / 830 ms

ECG Lead II

EGM1: Atp/Aring

ECG2: Vtp/Vring

Tests - Charge/Dump

Charge initiated...

1. Wait for charge to complete.

2. Retrieve Data.

Last Charge to Full Energy: 06-Jan-1994 04:42:11

Charge Time: 4.6 sec

Energy: 0.1 - 30 J

Last Capacitor Formation*: 06-Jan-1994 04:27:14

Buttons: RETRIEVE Data, ABORT Test, Print...