Kwiecień 2018 r.

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e-mail: serwis.medyczny@philips.com

Z poważaniem

Marta Łoniewska
Q&R Officer CEE North

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IntelliVue MX40

**Release B.05 or Later
Instructions for Use Errata**

This Errata is a supplement to the IntelliVue MX40 Instructions for Use Release B.05 and later. It contains important safety information. Attach this Errata to the first page of Chapter 6 to ensure that it is not misplaced and is stored with the Instructions for Use for ready reference.

Part Number 4535 647 99191
Printed in the U.S.A. April 2018
First Edition

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6 ECG and Arrhythmia Monitoring

ECG Safety Information

Warning The device provides QT and QTc interval change information; the clinical significance of the QT and QTc interval change information should be determined by a clinician. For more information, see the *QT Interval Monitoring Application Note*, p/n 452296278601.

For Paced Patients

- Warning**
- The output power of the MX40 and other sources of radio frequency energy, when used in the proximity of a pacemaker, can be sufficient to interfere with pacemaker performance. Due to the shielding effects of the body, internal pacemakers are somewhat less vulnerable than external pacemakers. However, caution should be exercised when monitoring any paced patient. In order to minimize the possibility of interference, position electrodes, electrode wires, and the MX40 as far away from the pacemaker as possible. Consult the pacemaker manufacturer for information on the RF susceptibility of their products and the use of their products with the MX40. See the *Patient Information Center Instructions for Use* for additional information on monitoring paced patients.
 - When an external pacemaker is being used on a patient, arrhythmia monitoring is severely compromised due to the high energy level in the pacer pulse. This may result in the arrhythmia algorithm's failure to detect pacemaker non-capture or asystole.
 - Pacemakers that create fusion beats (pace pulse on top of the QRS complex) cannot be detected by the monitor's QRS detector.
 - For paced patients who exhibit only intrinsic rhythm, the monitor can erroneously count pace pulses as QRS complexes when the algorithm first encounters them, resulting in missed detection of cardiac arrest. The risk of missing cardiac arrest can be reduced by monitoring these patients with the low heart rate limit at or slightly above the basic/demand pacemaker rate. A low heart rate alarm notifies you when the patient begins pacing. Proper detection and classification of the paced rhythm can then be determined.

Note — During defibrillation, monitoring may be temporarily interrupted or distorted. It may take several seconds for the ECG trace to reappear on the screen. After defibrillation, the device will continue to monitor as before; the device settings will not be affected.

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