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RegJoint™ is a tissue engineered implant indicated for arthroplasties in small joints in hands and feet. Specific target joints are metatarsophalangeal I-V, (MTP I-V), carpometacarpal I (CMC I), metacarpophalangeal II-V (MCP II-V) and proximal interphalangeal (PIP) joints.

RegJoint™ is manufactured of bioabsorbable poly-96L/4D-lactide copolymer fiber. It is a porous, disc like implant whose fibers offer an attachment base for the patient's fibroblasts and the porous structure enables the free flow of nutrients and cell metabolites in and out of the structure enhancing the cell welfare and allowing the neotissue ingrowth. It is essential to maintain implants porous structure by implantation. RegJoint™ sizers are available to help to choose the optimal size of RegJoint™ for each patient. RegJoint™ product and sizer selections are presented in Tables 1 and 2. RegJoint™ loses its initial strength during 15 to 24 weeks in vivo, with complete strength loss and resorption on average within 2-3 years, depending on the patient variables. RegJoint™ provides temporary support and guidance for the fibrotic tissue in-growth. It allows a gradual optimized replacement of the implant with connective tissue providing a flexible and durable pseudo joint with strong collagen framework.

Implanting RegJoint™ does not require any specific instrumentation. All used instruments need to be appropriately cleaned and sterilized before use. All sutures and other materials need to be sterile.


RegJoint™ is sterilized by gamma irradiation and is supplied sterile and ready for use. Reg-Joint™ shall not be re-sterilized by any method.

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Picture 6. Radial nerve and the deep branch of the radial artery.

Picture 7. Resection lines.

- nerve and the deep branch of the radial artery. Picture 6.
- The first carpometacarpal (CMC) joint is identified, the capsule released and opened dorsoradially
  - The first metacarpal is released carefully proximally to correct the prevailing adduction contracture
  - The resection of the proximal part of the first metacarpal is done by using an oscillating saw. The resection line of the proximal part of the first metacarpal is perpendicular with respect to the metacarpal. The resection of the cartilage surface of the trapezium is performed using a curette or an oscillating saw. Picture 7.
  - Other options are partial or complete trapezium resections without metacarpal resection. The extension of the resection is proportional to the joint laxity.
  - Perform synovectomy and revise osteophytes.
  - RegJoint™ sizers can be used to evaluate the correct size of the implant. While measuring the size do not pull the thumb. The ideal size completely covers the cortical edges but does not significantly exceed them (Picture 2). Make sure that the porous structure of the implant is retained after implantation.
  - RegJoint™ is fixed with transosseous resorbable sutures (e.g. PDS 2-0, with UCL needle ☺). Alternatively, the scaffold can be fixed into the joint capsule or with

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Adres: Scaffdex Oy