



Instructions for Use












---

MIAMI J<sup>®</sup> SELECT


		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	15
ES	Instrucciones para el uso	20
IT	Istruzioni per l'uso	25
NO	Bruksanvisning	30
DA	Brugsanvisning	34
SV	Bruksanvisning	38
EL	Οδηγίες Χρήσης	42
FI	Käyttöohjeet	47
NL	Gebruiksaanwijzing	51
PT	Instruções de Utilização	56
PL	Instrukcja użytkowania	61
CS	Návod k použití	66
TR	Kullanım Talimatları	70
RU	Инструкция по использованию	74
JA	取扱説明書	79
ZH	中文说明书	83
KO	사용 설명서	87

12

a


b


c


d


e


f

g

h

## DESCRIPTION OF SYMBOLS

Medical Device

Magnetic Resonance (MR) safe

### Product overview (Figure 1):

- a. Front
- b. Height Indicator Marks
- c. Chin Support
- d. Hook Landing Area
- e. Tracheal Opening
- f. Patient Compliance Lockout
- g. Height Adjustment Button
- h. Sternum Relief Knob
- i. Sternum Contact
- j. Back Panel
- k. Repeatable Fit Tabs (sold separately)
- l. Strap

## INTENDED PURPOSE

To provide gross immobilization to the cervical spine.

- X-ray and CT lucent

The device must be fitted and adjusted by a healthcare professional.

## INDICATIONS FOR USE

- Conditions requiring gross immobilization of the cervical spine. This may include:
  - C-Spine precaution for trauma patients
  - Immobilization for pre and post c-spine surgery

## CONTRAINDICATIONS

- Patients with a compromised airway or known spinal deformities such as ankylosing spondylitis.
- Patients with penetrating trauma injuries

## WARNING AND CAUTIONS

**Warning:** If an unstable fracture is suspected or unknown, with or without a sustained trauma, ensure additional spinal precautions are implemented to immobilize the spine.

**Warning:** Use of a cervical collar may increase intracranial pressure (ICP) through jugular venous compression.

**Warning:** Use of a cervical collar may increase complexity of airway management.

**Warning:** Cervical spine immobilization, including use of a cervical collar, has been associated with:

- Impaired respiratory effort and forced expiratory volume
- Pneumonia
- Aspiration
- Worsening of existing cervical spine injury
- Severe neurological deterioration in patients with ankylosing spondylitis
- Triggering of non-compliance or agitation

- Delay to definitive treatment
- Impaired physical examination/secondary survey

**Caution:** Cervical spine immobilization, including use of a cervical collar, has been associated with:

- Increased pain and discomfort, which may lead to increased spinal movement
- Impaired swallowing

**Caution:** Perfumes and harsh cleaning agents should not be used on or under the collar, as they may compromise the integrity of materials.

**Caution:** Use of a cervical collar for cervical spine immobilization is not recommended in situations where patients are awake, alert, not intoxicated, without neck tenderness or pain and without abnormal sensory or motor findings on examination.

**Caution:** Regular cleaning of the collar, pads, and the skin beneath along with inspection of the skin for any signs of irritation is required to reduce the risk of skin ulceration. Bedridden patients are at an increased risk for skin ulceration.

**Caution:** Federal law restricts this device to sale by, or on the order of, a physician.

#### **Precautions:**

- At least two persons are needed for initial collar placement: one to maintain the patient's head and neck in proper alignment, the second to fit the collar.
- Collar removal or any adjustments to the collar should be done with physician permission only and according to the physician's instructions.
- Unless otherwise specified by the physician, the patient should not remove the collar except to wash under it and change the pads.
- The patient will need a second person to help when cleaning the collar.

**These instructions do not supersede hospital protocol and/or direct orders of patient's physician.**

## **GENERAL SAFETY INSTRUCTIONS**

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

Patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions;
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

## **FITTING INSTRUCTIONS**

### **Sizing Tips**

- Miami J Select is an adjustable collar suitable for most adult and adolescent patients; some patient phenotypes (**Fig. 2**) and anatomical variations may not be accommodated.
  - The Miami J Select is suitable for Stout (**Fig. 2b**), Tall (**Fig. 2d**), Regular (**Fig. 2e**), and Short (**Fig. 2f**) patient phenotypes.

- The Miami J Select has not been tested in Super Short (**Fig. 2a**) and XSmall (**Fig. 2c**) patient phenotypes and should therefore not be used on these patients.
- Adjust to the height that fits most comfortably and maintains desired treatment position.
  - Height settings 2-5 will be used most frequently. Height settings 1, 6 and 7 will be required less frequently.
- Once patient can sit/stand following initial application, check spine is in desired treatment position and adjust height setting if necessary.
- To ensure proper sizing, it is a good idea to fit operative patients prior to surgery.
  - Collar size may be incorrectly identified in the presence of postoperative dressings and swelling.

### ***Application Instructions***

**Ensure that correct spinal protocol is followed throughout the application.**

1. Ensure the padding extends beyond the edge of the plastic and the Sternum Relief Knob is oriented vertically as shown in **Fig. 1**.
2. Slide the Back Panel behind the patient's neck and center it, noting arrow on panel pointing up (**Fig. 3**).  
**Note:** Long hair should be placed outside of the plastic.
3. Slide the Front of collar up the chest wall until the bottom of the tracheal opening is at the level of sternal notch. Chin Support may not contact chin at this point. (**Fig. 4**). Sides of the collar Front should be oriented up, off the trapezius, and toward the ears.
4. Holding the collar against the chest with the bottom of the tracheal opening at the level of the sternal notch, depress the Height Adjustment Button and manually adjust the Chin Support to achieve the desired height (**Fig. 5**).
5. When desired height is attained, release Height Adjustment Button to lock in place.
6. Verify that desired height is achieved on both sides of the brace. Depress Height Adjustment Button and make micro adjustments as necessary.
7. While holding the Front securely, place the sides of the collar Front within the sides of the Back Panel. Apply the Straps of the back panel to the Hook Landing Areas on the front. Tighten Straps alternately to an equal length on both sides. (**Fig. 6**)
8. If further height adjustment is required, depress Height Adjustment Button, and manually move to position. Use the Height Indicator Marks to record patient's collar height.
9. When desired position is attained, engage the Patient Compliance Lockout located behind the Height Adjustment Button by moving the lever to the left. (**Fig. 7**)
10. Straps must be aligned with the Hook Landing Areas. When patient is properly fit, there should be equal amounts of excess straps overhanging the Hook Landing Area. These may be trimmed (**Fig. 8**).

### ***Final Adjustments***

1. Optional Repeatable Fit Tabs can be applied to trimmed straps to provide a contact point for the patient to don and doff the collar to a repeatable position. To use, thread the strap end through the tab and attach to the hook molded into the tab (**Fig. 9**).

2. The blue Sternum Relief Knob adds a degree of freedom for the patient so that activities such as eating or talking can be performed with greater comfort. The knob turns 90 degrees to a horizontal position. **(Fig. 10)**
3. Upon physician recommendation, collar corners and edges may be modified or trimmed to relieve discomfort or pressure.

### **Final Fitting Checklist**

A properly applied Miami J Select will look like **Fig. 11**.

- Collar extends from mandible to just below sternal notch.
- Chin is centered comfortably in Chin Support. Chin should not extend over edge of Sorbatex™ pad, nor fall inside collar.
- Straps are aligned to the Hook landing area and are an equal length.
- Sides of Back should overlap sides of Front.
- Front of collar angled up toward ears. Lower plastic edge should not be resting on patient's clavicles nor digging into trapezius.
- No plastic touching skin. Sorbatex™ pads extend beyond all plastic edges.
- No slack or gaps in any of the straps.
- Tracheal opening and posterior vent are midline.
- Collar sits away from the neck. There should be a "finger" gap between the tracheal opening and the neck. If the collar is fit too closely/tightly to the neck, size up to the next taller size.

### **ACCESSORIES AND REPLACEMENT PARTS**

Please refer to the Össur Catalog for a list of available replacement parts or accessories.

### **CLEANING AND CARE (Fig. 12)**

- Remove the pads from the plastic shells.
- Hand wash the pads with mild soap and water. **DO NOT** use harsh detergents or bleach.
- Rinse. Wring out the excess water and squeeze in towel. Lay flat to air dry.
- Wipe the plastic collar shells clean with mild soap and water.
- Replace pads by attaching the grey side to the hook on the inside of the collar.

### **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

### **LIABILITY**

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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