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EN NeuroVac drainage system for closed wounds

NeuroVac flat silicone drainage, perforated

The product consists of a perforated, contrast-compatible flat drainage made of silicone with a round extension hose made of silicone.

NeuroVac^{Pro} flat silicone drainage, channeled

The product consists of a channeled, contrast-compatible flat drainage made of silicone with a round extension hose made of silicone and a trocar.

NeuroVac round silicone drainage, perforated


The product consists of a perforated, contrast-compatible round drainage made of silicone with a round extension hose made of silicone and an adapter.

NeuroVac^{Pro} round silicone drainage, channeled

The product consists of a channeled, contrast-compatible round drainage made of silicone with a round extension hose made of silicone, a trocar, and an adapter.

NeuroVac reservoir

The suction reservoir made of silicone is equipped with an anti-reflux valve.

General Information:

Before using this device read instructions carefully. The product may only be used according to the intended purpose as described in the following. The techniques described in these instructions are only a selection of the procedures which are medically possible. The physician having sole responsibility is free to perform the individual application steps according to a technique preferred by him. The manufacturer will not assume liability and/or warranty for damages caused by improper application and/or non-observance of the instructions.

Do not use a product when the individual package is broken or dirty or if the product itself looks abnormal in any way. Use the device prior to the Expiry Date on the package label. Do not resterilize.

NEUROMEDEX GMBH does not warrant for expired or resterilized products.

This medical product is intended for single use only. The reprocessing methods described in the specialist literature destroy the medical product or cause significant changes in the characteristics of the product and a safe application cannot be guaranteed.

Technical data

The components and technical data can be found on the label of the individual packaging. Further configurations can be found in the brochure.

Material

Flat / round drainage:	X-ray contrast-compatible silicone
Round hose:	silicone
Reservoir:	silicone / PP

Indications and fields of application

This drainage system is designed for being used as temporary low-vacuum wound drainage for closed wounds in the fields of abdominal surgery, plastic surgery, as well as after mamma reduction plastics.

Contraindications

Generally, incompatibility with one of the mentioned materials is regarded as contraindication.

Warnings

- Soft silicone materials are characterized by low tensile strengths. Thus, particular care has to be exercised when handling and positioning the product. Damages to the drainage hose due to scalpels or needles have to be avoided absolutely, as the same results in the tensile strength of the material being reduced significantly.
- If a drainage has been perforated accidentally, please do not use the same any longer.
- The maximum dwell time of this product is 6 days, as tissue ingrowth respectively agglutination may be the result.
- In order to be able to remove the drainage more easily at a later point in time, no part of the drainage hose – neither inside nor outside of the wounded area – must be rotated, pressed, overstretched, or fixed with suture instruments. **(This may result in material cracks and in the drainage being pulled off subsequently).**
- Drainages should be displaced or removed by exerting a uniform, light force and manually only. **Excessive pulling can result in the drainage being pulled off.**
- Drainage hoses and extensions must not be used with instruments of any type, as the same could affect the material (cracks, deformations), which in turn could result in the drainage being pulled off.
- In order to prevent tissue ingrowth and infections, the drainage should be removed as soon as wound secretion is no longer removed.
- **The material being pulled of when removing the drainage can also occur in case of tissue ingrowth around drainage and perforations, if the dwell time has been extended accordingly.** Thus, the physician should monitor the development of the wound healing process with the patient carefully. If tissue ingrowth has been tolerated, removing the drainage surgically may be required.
- Drainages, trocars, connectors, and reservoirs are designed exclusively for being used with an individual patient and should be disposed of after being used.
- In any case, you have to prevent the drainage from becoming agglutinated or clogged with blood clots, or the reservoir from being filled completely with secretion, as the wound would no longer be drained this way. If the drainage clogs, the same can be sucked free by connecting a suction pump to the discharge aperture of the reservoir or directly to the drainage. If the reservoir is not emptied before it is filled completely, the likelihood of recontamination via the anti-reflux valve increases.
- Before starting the operation the patient should be informed on risks and complications that may result from a possible rejection reaction of his body - caused by material intolerances.

Handling and application

Please observe the following directives on positioning the drainage:

In order to prevent the silicone drainage from being contaminated with particles, the drainage should be taken from its packaging immediately before being used only. In doing so, please check the packaging for possible damages.

Carefully check the drainage and do not use the same if you detect any damages to the product or the packaging.

Drainage

1. The drainage is inserted into the affected wound area, drained, and removed by means of point incision through the adjacent skin area. The drainage may be cut to the required dimensions.
2. Please ensure that the perforated / channeled section of the drainage is located completely inside the wound area. The perforations / channeling should not be within the subcutaneous tissue (risk of subcutaneous hematomas). Insert the round hose through the incision spot until the non-perforated / non-channeled drainage area becomes visible at the skin surface.
3. All drainage apertures have to be inside the wound area. In order to obtain the vacuum, the wound has to be closed in an air-tight manner. If air enters at any position, this has to be corrected immediately, as the system will not work properly otherwise.
4. The drainage can be fixed with suture material slings. **The drainage must not be perforated with a needle and the suture material must not be guided through the drainage in any case.**

Reservoir

1. Immediately before being used the reservoir should be checked for patency (see item 3). If the reservoir cannot be filled with air completely or of this procedure takes more than 5 seconds, please do not use this reservoir.
2. For simple connection the drainage hose respectively the adapter will be pushed over the access aperture #1 of the reservoir as far as it will go.

- In order to generate the desired vacuum inside the reservoir, activate the discharge aperture and press-out as much air as possible from the reservoir. The completely drained reservoir is now sealed by inserting the sealing cap into the discharge aperture as far as it will go.
- In order to drain the secretion, the sealing cap is pushed from the discharge aperture and the reservoir is positioned above a collection basin in a slightly tilted position. By pressing together the reservoir the secretion is pressed out now. When the reservoir is drained completely, the vacuum inside the reservoir has to be re-generated before the same can be reused (see item 3).
- In order to generate constant, uniform suction, the reservoir should be positioned on patient level or below the drainage outlet. The clip is designed to attach the reservoir.
- The drainage can also be connected to an external suction pump or a wall-mounted suction unit. For this, using an extension hose is required, however.

Disposal

After use discard the product according to the regional or national waste regulations for infectious materials.

Exclusion of liability

We declare on our own responsibility that the product is in compliance with the applicable provisions of the directive 93/42/EEC for medical devices. In consideration of the biological differences of patients, the efficiency of the product cannot absolutely be guaranteed.

As the conditions of use, the establishment of diagnosis and indication as well as the application and use of the product are beyond our control, we cannot guarantee the success of the product or exclude the possibility of side effects.

These instructions have been written in German originally and have been translated into English. In case of discrepancies in interpretation the German version has priority.

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