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8840 N'Vision Programmer Printouts

[Print Report]

```

Print Report
8840 SynchroMed II B
Programming Session

Current Pump Settings
Examined At
 09/18/2008  MMDD/YYYY  10:39
Last Change
 08/18/2008  MMDD/YYYY  12:11

Pump in Safe State
Reset Occured
Patient Information
          
```

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[Event Log]


```

Event Log
System
MM/DD/YYYY  h:m
8/18/2008  12:13  (07)
Pump in Safe State
8/18/2008  12:13  (28)
Reset Occured - Low Battery
8/18/2008  12:13  (01)
Reset Occured
          
```

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
8840 N'Vision Programator wydruków

[Drukuj raport]

```

Print Report

8840 SynchroMed II B
Programming Session

Current Pump Settings
Examined At      01/14/2009  MM/DD/YYYY  13:39
Last Change     01/14/2009  MM/DD/YYYY  13:39
ERI Occurred
Schedule to Replace the pump by
04/10/2009  MM/DD/YYYY

Patient Information
First Name
XXXXXXXXXX
    
```

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```

Infusion Mode
Simple Continuous
Drug
Dose Per Day
1 LIORESAL
1,212.1 ug/day

ERI Occurred
01/10/2009  MM/DD/YYYY
    
```

[Dziennik zdarzeń]

```

Event Log
System
MM/DD/YYYY  h:m
01/10/2009  20:52  (OD)
ERI Occurred
    
```

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Medtronic

**INDICATIONS, DRUG STABILITY,
AND EMERGENCY PROCEDURES**

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Drug stability 5

Emergency procedures 7

Morphine intrathecal/epidural overdose 7

Baclofen injection overdose 9

Baclofen injection underdose/withdrawal 11

Emergency procedure to empty the pump reservoir 12

**Refer to the appropriate information for prescribers booklet for contraindications, warnings, precautions, adverse events summary, individualization of treatment, patient selection, use in specific populations, and component disposal.
Refer to the device implant manual for device description, package contents, device specifications, and instructions for use.**

A12510_001TOC.fm 6/14/06 3:36 pm
UC200xxxxxx EN
4.5 x 6 inches (114 mm x 152 mm)

Medtronic Confidential
NeuroIFP_R00

4 English Indications, drug stability, and emergency procedures 2006-06

MA12510A005 Rev A

Indications

The SynchroMed II, SynchroMed EL, and SynchroMed Infusion Systems and the IsoMed Constant-flow Infusion System are indicated for use when patient therapy requires the chronic infusion of drugs referred to in Table 1.

Drug stability

Physicians prescribing the SynchroMed II, SynchroMed, or SynchroMed EL Infusion Systems or the IsoMed Constant-Flow Infusion System for use with the drugs listed in Table 1 must be familiar with the indications, contraindications, warnings, precautions, adverse events, dosage and administration information, and screening procedures described in the drug labeling. Each system includes (at a minimum) a pump and a catheter. Testing has indicated that the drugs in Table 1 are stable and compatible with the infusion systems listed in the table. Refer to the appropriate drug labeling for complete prescribing information, including indications, contraindications, warnings, precautions, and adverse events.

**Table 1. Stability^a of drugs approved for use with
Medtronic implantable infusion systems**

	SynchroMed II	SynchroMed and SynchroMed EL	IsoMed
Floxuridine (FUDR) (20 mg/mL)	56 days	28 days	27 days
Baclofen injection (0.5 mg/mL) (2 mg/mL)	180 days	90 days	35 days
Methotrexate (5 mg/mL)	56 days	28 days	28 days
Morphine sulfate sterile solution (preservative-free) (25 mg/mL)	180 days	90 days	90 days
Morphine hydrochloride sterile solution (preservative-free) (25 mg/mL)	180 days	90 days	90 days

**Table 1. Stability^a of drugs approved for use with
Medtronic implantable infusion systems (continued)**

	SynchroMed II		SynchroMed and SynchroMed EL		IsoMed
Ziconotide sterile solution (preservative-free)	Initial fill^b	Refill	Initial fill	Refill	—
• 25 µg/mL, undiluted	14 days	60 days	14 days	60 days	
• 100 µg/mL, undiluted	NA	60 days	NA	60 days	
• 100 µg/mL, diluted	NA	40 days	NA	40 days	

^a Stability is defined as 90% of initial concentration.

^b For a pump that has not been previously filled with preservative-free ziconotide sterile solution, only the undiluted 25 µg/mL formulation of preservative-free ziconotide sterile solution can be used for the initial pump fill. Refill of the pump with this drug should be completed within 14 days to avoid underdosing the patient. At initial fill of a new pump, some of the drug is lost due to two factors that do not occur upon subsequent refills: adsorption on internal device surfaces, such as the titanium, and by dilution in the residual space of the device. Refer to the drug labeling for additional information.

— Drug is not approved for this pump.

Emergency procedures

Morphine intrathecal/epidural overdose

Consult the patient's medical record or with the patient's physician to confirm the drug or drug concentration within the pump reservoir.

Symptoms

Respiratory depression with or without concomitant central nervous system depression (ie, dizziness, sedation, euphoria, anxiety, seizures, respiratory arrest).

Actions


Figure 1. Morphine intrathecal/epidural overdose emergency procedures.

- 1 Preservative-free morphine sulfate sterile solution manufacturer's package insert.
- 2 Naloxone hydrochloride manufacturer's package insert.
- 3 Refer to the drug manufacturer's package insert for a complete list of indications, contraindications, warnings, precautions, adverse events, and dosage and administration information.
- 4 Use a 25-gauge needle for withdrawal from a SynchroMed or SynchroMed EL catheter access port. Use a 24- or 25-gauge needle for withdrawal from a SynchroMed II or IsoMed catheter access port.

Baclofen injection overdose

Consult the patient's medical record or with the patient's physician to confirm the drug or drug concentration within the pump reservoir.

Symptoms

Drowsiness, lightheadedness, dizziness, somnolence, respiratory depression, seizures, rostral progression of hypotonia, and loss of consciousness progressing to coma.

There is no specific antidote for treating overdoses of baclofen injection. However, anecdotal reports suggest that intravenous physostigmine may reverse central side effects, notably drowsiness and respiratory depression.

Actions


Figure 2. Baclofen injection overdose emergency procedures.

- 1 Refer to the drug manufacturer's package insert for a complete list of indications, contraindications, warnings, precautions, adverse events, and dosage and administration information.
- 2 Müller-Schwefe G, Penn RD. Physostigmine in the treatment of intrathecal baclofen overdose: report of three cases. *J Neurosurg.* August 1989;71:273-275.
- 3 Use a 25-gauge needle for withdrawal from a SynchroMed or SynchroMed EL catheter access port. Use a 24- or 25-gauge needle for withdrawal from a SynchroMed II or IsoMed catheter access port.

Baclofen injection underdose/withdrawal

Consult the patient's medical record or with the patient's physician to confirm the drug or drug concentration within the pump reservoir.

Symptoms of underdose

Pruritus without rash, hypotension, paresthesia, fever, and altered mental state.

Symptoms of withdrawal

Exaggerated rebound spasticity and muscle rigidity, rhabdomyolysis, and multiple organ failure. The condition may resemble autonomic dysreflexia, sepsis, malignant hyperthermia, and neuroleptic-malignant syndrome.

Actions


Figure 3. Baclofen injection underdose/withdrawal emergency procedures.

1 Refer to the drug manufacturer's package insert for a complete list of indications, contraindications, warnings, precautions, adverse events, and dosage and administration information.

2 Coffey RJ, Edgar TS, Francisco GE, et al. Abrupt withdrawal from intrathecal baclofen: recognition and management of a potentially life-threatening syndrome. *Arch Phys Med Rehabil.* 2002;83:735-741.

Emergency procedure to empty the pump reservoir

Equipment

- 22-gauge needle
- 20-mL syringe(s)
- 3-way stopcock
- Antiseptic agent

1. Locate the pump by palpation. The reservoir fill port is located in the CENTER of the pump.
2. Prepare the injection site by cleansing the area using antiseptic agent.
3. Gently insert the 22-gauge needle into the center of the reservoir fill port septum until the needle touches the needle stop. If the clinician encounters resistance during needle insertion, the clinician should reassess placement. Do not force the needle. Hitting metal or the feel of abnormal resistance during the procedure may be an indication that the needle is not in the center of the reservoir fill port septum.

△ Cautions:

- The IsoMed pump reservoir contents are under significant pressure. To prevent the reservoir contents from being ejected, do not use an open syringe when emptying the pump.
 - When removing a vesicant or cytotoxic drug, care must be taken to prevent spillage or leakage of the drug into adjacent tissue.
4. If desired, use a 3-way stopcock. Withdraw fluid from the pump reservoir using passive backflow and gentle negative pressure. Empty the pump reservoir until back-flow has stopped and wait five seconds to ensure all fluid is removed and the pump reservoir is empty. Depending on pump reservoir volume, more than one syringe may be needed to empty the pump.
 5. Remove the needle from the reservoir fill port.
 6. Record in the patient chart the amount of fluid emptied from the pump reservoir.

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Medtronic®, IsoMed®, and SynchroMed® are registered trademarks of Medtronic, Inc.
All Rights Reserved
Medtronic International Ltd.
© Medtronic, Inc. 2006