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**Grzegorz Cessak**

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# Report Form

## Field Safety Corrective Action

Medical Devices Vigilance System  
(MEDDEV 2.12/1 rev 5)

1. Administrative information	
<b>Destination</b>	
Name of national competent authority (NCA) MHRA	
Address of national competent authority 151 Buckingham Palace Road, Victoria, London, SW1W 9SZ	
Date of this report 06.02.2012	
Reference number assigned by the manufacturer 06.02.2012	
Incidence reference number and name of the co-ordinating national competent authority (if applicable)	
Identify to what other national competent authorities this report was also sent EU: Austria, Belgium, Bulgaria, Switzerland, Czech Republic, Germany, Denmark, Estonia, Spain, Finland, France, Great Britain, Greece, Hungary, Ireland, Italy, Lithuania, Luxembourg, Latvia, Malta, Netherlands, Poland, Romania, Sweden, Slovenia, Slovakia, Croatia, Turkey. Non-EU: Argentina, Bangladesh, Belarus, Brazil, Bosnia, China, Colombia, Costa Rica, Egypt, Guatemala, Hong Kong, Israel, Kenya, Korea, Laos, Macedonia, Malaysia, Mexico, Nicaragua, Pakistan, Panama, Peru, Philippines, Serbia, Singapore, Ukraine, Uruguay, Venezuela, Vietnam.	
2 Information on submitter of the report	
Status of submitter	
<input checked="" type="checkbox"/> Manufacturer	
<input type="checkbox"/> Authorised representative within EEA	
<input type="checkbox"/> Others (identify the role):	
3 Manufacturer information	
Manufacturer name Alcon Laboratories Inc	
Manufacturer's contact person Stephen Beekman, Manager Post Market Compliance	
Address 6201 South Freeway	
Postal code 76134-2099	City Fort Worth, Texas
Phone +00 1817 615 5491	Fax +00 1817 916 9087
E-mail stephen.beekman@alconlabs.com	Country USA
4 Authorised representative information	
Name of the authorised representative Alcon Laboratories (U.K.) Ltd.	
The authorised representative's contact person Hilde Viroux	
Address Pentagon Park, Boundary Way, Hemel Hempstead	
Postal code HP2 7UD	City Hertfordshire
Phone +44 1442 341 175	Fax +44 1442 341 280
E-mail authorised.representative@alconlabs.com	Country United Kingdom
5 National contact point information	
National contact point name Alcon Polska Sp. z o.o.	

Name of the contact person Joanna Morgan	
Address Żelazna 28/30	
Postal code 00-832	City Warsaw
Phone 0-22/ 820 34 50	Fax 0-22/ 820 34 56
E-mail joanna.morgan@alconlabs.com	Country Poland

### 6 Medical device information

Class	
<input type="checkbox"/>	AIMD Active implants
<input type="checkbox"/>	MDD Class III
<input checked="" type="checkbox"/>	MDD Class IIb
<input type="checkbox"/>	MDD Class IIa
<input type="checkbox"/>	MDD Class I
<input type="checkbox"/>	IVD Annex II List A
<input type="checkbox"/>	IVD Annex II List B
<input type="checkbox"/>	IVD Devices for self-testing
<input type="checkbox"/>	IVD General
Nomenclature system (preferable GMDN) GMDN	
Nomenclature code 47410	
Nomenclature text Anterior-chamber Intraocular Lens, Phakic	
Commercial name/brand name/make ACRYSOF® CACHET™ Phakic Lens	
Model number L12500, L13000, L13500, and L14000	
Serial number(s) and/or lot/batch number(s) All	
Software version number (if applicable) n/a	
Manufacturing date/expiry date (if applicable) All	
Accessories/associated device (if applicable) n/a	
Notified body (NB) ID- number 0123	

### 7 Description of FSCA

Background information and reason for the FSCA  
The ACRYSOF® CACHET™ Phakic Lens is a UV-absorbing anterior chamber meniscus lens designed to provide refractive correction in the treatment of moderate to high myopia. The CACHET® Lens is a CE marked medical device marketed by Alcon since 2008 in approximately 70 countries. Recently, data from four and five year clinical studies conducted by Alcon has been analyzed. The additional information gained through the studies allows an update to the Directions for Use (DFU). These updates are primarily associated with:

- Clinical data related to subject accountability
- Mean Annualized % Central Cell Density Loss
- Age range brackets are being split and minimum endothelial cell density (ECD) requirements are increasing for specific age groups
- Warnings will be updated to describe guidance for establishing a baseline to monitor chronic Endothelial Cell Loss (ECL), recommendations for when to conduct more frequent monitoring, and guidance regarding when to consider IOL explantation

ECL is the reduction in the number of cells of the corneal endothelium and is a known risk associated with the placement of a phakic lens. There are no clinical symptoms of ECL until the ECL becomes severe enough to degrade optical performance and can compromise visual acuity. While information to mitigate the risk of ECL is already included in the CACHET® DFU, the additional clinical study data allows an update to the DFU to provide additional information and guidance in order to minimize the risks associated with ECL.

This FSCA is not related to product quality issues or adverse events reported from marketed product.

Description and justification of the action (corrective/preventive)

To communicate these updates to customers, Alcon feels that a customer communication is necessary. Because the updates to the DFU are safety related and the device is CE marked, a customer letter meets the definition of a Field Safety Notice (FSN) per MedDev 2.12-1 Rev 6 (December 2009) and a Field Safety Corrective Action is required.

This FSCA is not related to product quality issues or adverse events reported from marketed product.

Advice on actions to be taken by the distributor and the user

Physicians should become familiar with the updates to the DFU, specifically those updates associated with ECD and ECL.

Attached please find

- Field Safety Notice (FSN) in English  
 FSN in national language  
 Others (please specify):

Time schedule for the implementation of the different actions

FSNs will be distributed within countries where the ACRYSOFF® CACHET™ Phakic Lens is marketed based upon regulatory approval for the DFU update and/or the regulatory approval of the FSN.

These countries within the EEA and Switzerland are affected by this FSCA

Within EEA and Switzerland:

- AT    BE    BG    CH    CY    CZ    DE    DK    EE    ES  
 FI    FR    GB    GR    HU    IE    IS    IT    LI    LT  
 LU    LV    MT    NL    NO    PL    PT    RO    SE    SI  
 SK

Candidate Countries:

- HR    TR

All EEA, Candidate Countries and Switzerland

Others:

These countries outside the EEA and Switzerland are affected by this FSCA

Argentina, Bangladesh, Belarus, Brazil, Bosnia, China, Colombia, Costa Rica, Egypt, Guatemala, Hong Kong, Israel, Kenya, Korea, Laos, Macedonia, Malaysia, Mexico, Nicaragua, Pakistan, Panama, Peru, Philippines, Serbia, Singapore, Ukraine, Uruguay, Venezuela, Vietnam

#### 8 Comments

I affirm that the information given above is correct to the best of my knowledge.

  
Signature

Joanna Morgan  
Name

Warsaw  
City

06.02.2012  
Date

*Submission of this report does not, in itself, represent a conclusion by the manufacturer and/or authorized representative or the national competent authority that the content of this report is complete or accurate, that the medical device(s) listed failed in any manner and/or that the medical device(s) caused or contributed to the alleged death or deterioration in the state of the health of any person.*