zimmer

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ZMR[®] Porous Revision Hip Prosthesis

Surgical Technique

ZMR Porous Revision Hip Prosthesis Surgical Technique

Table of Contents


Preoperative Planning	2
Objectives Of Preoperative Planning	2
Assessment of Bone Loss and Component Selection	3
Determination of Leg Length	3
Determination of Femoral Offset and Abductor Muscle Tension	3
Templating For The ZMR Porous Revision Hip System	4
Determining The Type, Size, Build-up, And Offset Of The Modular Proximal Body	5
Determining The Modular Distal Stem Size And Length	6
Surgical Technique	8
Distal Canal Preparation	8
Proximal Femur Preparation	10
Trial Reduction	14
Implantation	17
Femoral Head Assembly	20
Wound Closure	20
Postoperative Management	21
Implant Removal	21
What ZMR Instruments do you need	24
Instruments	25
Instrumentation Options	32

Objectives of Preoperative Planning

The objectives of preoperative planning include:

1. Assess amount of bone loss to:
 - Select type(s) of implant(s) appropriate to reconstruct the femur.
2. Ensure that the proximal body will be fully supported in the area of the plasma spray.
3. Obtain the anticipated component size necessary to provide structural stability.
4. Establish parameters of joint kinematic restoration that include:
 - Determination of leg length.
 - Restoration of offset for abductor muscle tensioning.
5. Assess the acetabulum to determine if any acetabular reconstruction is needed, and consider the potential impact of the acetabular reconstruction on the femoral side.

WARNING: The standard junction **ZMR Hip System** should only be used when full proximal support will be achieved in the area of the plasma spray. This is necessary because without full proximal support, the mid-stem junction is vulnerable to fracture. If such proximal support cannot be achieved, evaluate the use of the **ZMR XL**.

Note: The Porous Coated areas shown must be in contact with the metaphyseal femur.

Assessment of Bone Loss and Component Selection

To select the appropriate type of implant in femoral revision surgery, a number of factors must be evaluated. Reconstructing the femur based on the amount of femoral bone loss and the status of the remaining bone is important in determining the appropriate prosthesis.

The ZMR Porous Revision Hip System is a versatile system that can address a wide range of femoral revision needs. It offers the surgeon multiple fixation options. These options include proximal fixation, and combined proximal and distal fixation (extensive fixation). The implant choice will depend on the individual patient and the surgeon's fixation preference.

Determination of Leg Length

Preoperative determination of leg length is essential for the restoration of the appropriate leg length during surgery. In femoral revision surgery, correction of leg length discrepancy is usually necessary because of bone and soft tissue changes resulting from the failed prosthesis.

An anterior/posterior (A/P) pelvic radiograph often provides enough detail of leg length inequality to proceed with surgery. If more information is needed, a full-length femoral x-ray view may be helpful. From the clinical examination and radiographic information on leg length, the needed correction, if any, can be determined.

If leg length is to be maintained or minimally increased, it is usually possible to perform the operation successfully without osteotomy of the greater trochanter, unless extensive exposure of the acetabulum or femur is necessary. However, if there is some major anatomic abnormality, osteotomy of the greater trochanter may be helpful and safer.

In the unusual situation where the limb is to be significantly shortened, subtrochanteric osteotomy, or osteotomy and advancement of the greater trochanter are mandatory. If the limb is shortened without one of these techniques, the abductors will be lax postoperatively, and the risk of dislocation will be high. Also, gait will be compromised by the laxity of the abductors.

Determination of Femoral Offset and Abductor Muscle Tension

After establishing the desired leg length requirements, abductor muscle tensioning through femoral offset must be considered. Restoring adductor tension places the abductor muscles at their optimal mechanical advantage, and results in improved functioning with less chance of dislocation and limp.

When the patient has a very large offset between the femoral head center of rotation and the line that bisects the medullary canal, the insertion of a femoral component with a lesser offset will, in effect, medialize the femoral shaft. To the extent that this occurs, laxity in the abductors will result.

Using the ZMR Porous Revision Templates on the preoperative x-ray film allows the surgeon to estimate the amount of offset needed for the new femoral component to either improve or reproduce the patient's anatomical offset.

WARNING: The standard junction ZMR Hip System should only be used when full proximal support will be achieved in the area of the plasma spray. This is necessary because without full proximal support, the mid-stem junction is vulnerable to fracture. If such proximal support cannot be achieved, evaluate the use of the ZMR XL.

Templating for the ZMR Porous Revision Hip System

Preoperative planning for the insertion of the ZMR Porous Revision Femoral Component requires at least two views of the involved femur: an A/P view of the of the pelvis centered on the pubic symphysis, and a Lowenstein lateral view on an 11" x 17" cassette. Both views should show the full length of the femur. In addition, it may be helpful to obtain an A/P view of the involved side with the femur internally rotated. This compensates for naturally occurring femoral version, and provides a more accurate representation of the true medial-to-lateral dimension of the metaphysis.

When templating, magnification of the femur will vary depending on the distance from the x-ray source to the film, and the distance from the patient to the film. The ZMR Porous Revision Templates use 20 percent magnification, which is near the average magnification on most clinical x-rays. Large patients and obese patients may have magnification greater than 20 percent because osseous structures are farther away from the surface of the film. Likewise, smaller patients may have magnification less than 20 percent. If necessary, to better determine the magnification of any x-ray film, use a standardized marker at the level of the femur. (Templates of other magnification can be obtained as a special order by contacting your Zimmer Sales Representative.)

Preoperative planning is important in choosing the acetabular component if it is also being revised. Acetabular preoperative planning is beyond the scope of this document, but must be considered in conjunction with the use of the ZMR Porous Revision Femoral Prosthesis.

The objectives in templating the ZMR Porous Revision Femoral Component include:

1. Determining the type, size, position, length, and offset of the modular proximal body.
2. Choosing the type, size, length, and position of the distal stem.
3. Developing a plan for adjunctive support of the proximal body superior to the mid-stem modular junction.

The ZMR Porous Revision Templates include separate proximal body and distal stem templates. These templates are used together to help determine the final implant component.

To fit the Spout, Cone, and Calcar Bodies, separate templates are available for each of these ZMR body styles.

The Spout Bodies are designed with a medial spout (curve) to provide fill for proximal fixation.

The Cone Bodies, which are conical in shape, are designed for accurate fit at the metaphyseal/ diaphyseal junction. They provide a wide range of version adjustments.

The Calcar Bodies are conical in shape. They have a collar that is designed to rest on the medial area of the femur for accurate fit at the metaphyseal/ diaphyseal junction.

The build-up heights available with these proximal bodies, combined with femoral head-neck length options help restore leg length. In addition, the body offset options, in combination with the modular femoral heads, assist in adjusting lateral offset and abductor muscle tensioning.

Another way to restore or compensate for offset is osteotomy and advancement of the greater trochanter to increase tension in the abductor muscles. In addition, the use of special acetabular liners can change the placement of the joint center of rotation.

The distal stem templates for the ZMR Porous Revision Hip Components are used for splined and fully porous coated options. The polished, splined distal stems provide distal augmentation to control rotational forces. The porous-coated distal stems provide distal fill for axial and rotational control with porous coating for implant/bone fixation.

The distal stem templates show the modular junction level as well as the various stem lengths and diameters. The templates also indicate those stem sizes that have both straight and bowed options. The diameter measurements on the splined stem templates include the height of the splines. The diameter measurements of the porous stem templates are the approximate measurements taken over the porous coating.

Determining the Type, Size, Build-up, and Offset of the Modular Proximal Body

Once the hip joint center of rotation has been established by templating the acetabulum, the modular proximal porous femoral body should be templated first. This allows the selected proximal body to determine where the mid-stem modular junction level will be prior to templating the distal stem component.

Selection of the proximal body type is determined in part by the femoral bone deficiency as assessed on the A/P radiograph. When the femoral tube is intact, the Spout Body can be used if proximal fill is desired. However, the Cone Body can also be used, particularly if special consideration must be given to version. When there is medial segmental femoral bone loss, either the Calcar Body or the Cone Body can be used. The various body heights allow the prosthesis to be deeply sealed in the canal so it can be anchored in viable bone stock.

Each porous body is available in sizes that increase progressively in the A/P and M/L dimensions. Both A/P and lateral radiographs are essential in determining the appropriate sized implant. Lay the selected proximal body template over the A/P x-ray film to determine the fit (Fig. 1). At all times consider the position of the femoral head center of rotation relative to the hip center in assessing leg length and offset rotation. Tight apposition of the implant to bone at the metaphyseal/diaphyseal junction is important with all the proximal bodies; however, it is especially crucial with the Cone and Calcar Bodies because they rely on the flare of the distal portion of the proximal body to provide stability against axial load. Use the lateral radiograph to further assess the fit of the proximal body (Fig. 2).


Fig. 1


Fig. 2

The A/P x-ray view is useful in determining the correct leg length and offset needed for the reconstruction. The various build-up options (35mm, 45mm, 55mm) available with the modular bodies and femoral heads will allow for adjustment of leg length. The Proximal Bodies are available with reduced (36mm), standard (40mm) and extended (46mm) offsets. The reduced offset bodies offer medial translation at 4mm. The extended offset bodies offer lateral translation of 6mm. This allows for an offset increase of 6mm without changing the vertical height or leg length. The femoral head lengths will also affect leg length and offset (Table 1).

Table 1: Adjustments available with the implant for restoring joint kinematics

Adjustments	Body	Head
Offset	Using reduced, standard and extended offset bodies will change offset only (up to 10mm).	Using varying head lengths will change both offset and leg length (up to 10mm in 2.5mm increments)†
Leg Length	Using 35mm, 45mm, and 55mm body heights will change leg length only (up to 20mm in 10mm increments).	Using varying head lengths will change both leg length and offset (up to 10mm in 2.5mm increments).
Version	Adjustment is possible by rotating body on the stem for Calcar and Cone Body options.	No adjustment is possible.

† Using VerSys® Hip System 12/14 cobalt chrome heads

Determining the Modular Distal Stem Size and Length

Use the distal stem templates in conjunction with the selected proximal body template to estimate the appropriate size and length of the distal stem. With the proximal body template still in position on the A/P x-ray film, or marked through the template holes, superimpose the distal stem template on the isthmus, making sure it aligns with the +0mm red line on the proximal body template. This will correctly align the mid-stem modular junction of the proximal body and distal stem.

Estimate the size and length of the distal stem (Fig. 3). The distal stem should fill, or nearly fill, the medullary canal in the isthmus area on the A/P x-ray view. To provide for adequate stability, the stem length should be sufficient to engage at least 50mm of intact diaphyseal bone.

Next, overlay the proximal body and distal stem templates on the lateral radiograph to assess fill and appropriate stem length while taking the anterior bow of the femur into consideration (Fig. 4). Again, it is important that the distal stem template be aligned correctly with the +0mm red line on the proximal body template. In a femur with a pronounced anterior bow, the bowed stem should be considered. Another option is to use a shorter straight stem if there is adequate femoral canal engagement. The Spline Stems have a coronal slot, which facilitates insertion of a longer straight stem.

It is also important to have adequate stem length when bypassing any femoral defects. When bypassing a cortical defect, the implant should extend past the defect by a minimum of two and one-half times the measured canal diameter to provide adequate support and to decrease the risk of femoral fracture.


Fig. 3


Fig. 4

Surgical Technique Steps

1	Distal Canal Preparation	8
	A Straight Stems	9
	B Bowed Stems.	10
2	Proximal Femur Preparation	10
3	Trial Reduction	14
4	Implantation.	17
	A Back Table Assembly	17
	B Loose Assembly Technique	20

1 Distal Canal Preparation

The distal femoral canal must be prepared to accept the distal femoral stem. The Spline Stems, when used in conjunction with the Spout Body, provide for distal adjunct stability, particularly rotational stability. This facilitates the proximal fixation of the Spout Body. The Porous Stems provide for distal fixation that can be used with any of the proximal bodies for an extensive fixation option.

A – Straight Stem Femoral Canal Preparation

When preparing the distal femoral canal for a straight distal stem component, use the straight intramedullary reamers (Fig. 5). Another option is to use the T-handle with the straight reamers for hand reaming. It is important to ream to the proper depth to assure adequate seating of the implant. The reaming depth should be greater than the preoperatively planned length. This will accommodate any change to the planned components.

Begin reaming with a straight reamer that is 4mm or 5mm smaller than the anticipated prosthesis size. Sequentially increase the reamer size by 0.5mm increments, making sure that each reamer is fully advanced to its appropriate depth and centered in the medullary canal. The appropriate reamer depth should be determined preoperatively when assessing the length of implant required. However, to accommodate any possible changes to the planned components, the reaming depth should be greater than the preoperatively planned length. The VerSys® IM Reamer depth markings **DO NOT** directly correspond to the ZMR stem length options. These markings are provided as guides to be referenced against anatomical land marks in the proximal femur (e.g. tip of greater trochanter) to aid in achieving proper depth of reaming. Avoid varus positioning of the reamers. If the greater trochanter or lateral neck tends to push the reamer medially, use a small conical reamer (trochanteric router) to clear the obstruction and then resume the straight distal reaming. Ream until adequate cortical contact has been achieved. When using the Spline Stem option, underreaming of the canal by 0.5mm less than the implant (e.g., ream to 14.5mm to implant a size 15mm distal stem) will provide for additional apposition of the

distal splines with the femoral canal. The Spline Stem can be reamed line-to-line based on surgeon preference.

If a straight porous stem is chosen to provide axial stability, the canal should be underreamed by 0.5mm or reamed line-to-line. It is good practice to measure the diameter of the last reamer and compare it to the diameter of the actual prosthesis. Use an accurately calibrated caliper to perform this task.

Ultimately, the reaming decision is based on the patient's anatomy, length of stem, bone quality, and the surgeon's judgment and experience.

Fig. 5

B – Bowed Stem Femoral Canal Preparation

When preparing the femoral canal to accept a bowed stem, flexible reamers must be used because they more closely follow the natural bow of the femur. The *Pressure Sentinel*® Intramedullary Reaming System has flexible reamers that result in lower intramedullary pressure levels and lower temperature levels that are at or below those of competitive systems.¹

Caution: Tactile feedback with these flexible reamers is different than straight reamers. Fluoroscopy can be used, if desired, to monitor the reaming process.

Insert a ball-tipped guide wire into the canal, and begin reaming 2mm or 3mm below the anticipated implant size (Fig. 6). Sequentially increase the reamer size by 0.5mm increments, making sure the reamer is fully advanced by passing through the femoral diaphysis. Ream until the desired canal diameter has been achieved. When reaming for the bowed Spline Stem, line-to-line to 1.0mm overreaming techniques can be used, depending on the patient's anatomy, bone quality, and surgeon judgment. The Spline Stem has a distal coronal slot design that makes the distal stem more flexible to accommodate the bow of the femur. When reaming for the bowed Porous Stem, it may be necessary to overream by as much as 1.0mm to accommodate the bow of the femur. In revision surgery, where the primary axial stability is achieved distally with an extensively porous coated stem, it may be necessary to ream line-to-line or overream 0.5mm to achieve this goal. The bowed Porous Stems do not have a coronal slot and are not as accommodating to the bow of the femur as slotted stems. Therefore, overreaming will facilitate the insertion of these stems.

Fig. 6

2 Proximal Femur Preparation

Bone must be cleared from the proximal femur to allow for intimate apposition of the proximal body with bone. Regardless of which body style is used, proper preparation of the body conical region helps secure the body at the femoral metaphyseal/diaphyseal junction and provides resistance to axial load. When using the Spout Body, additional bone must be removed from the medial calcar region of the femur to match the spout geometry. The Spout Body is designed for use when medial fill in the metaphysis is desired.

When preparing for a Calcar Body, use the Osteotomy Guide as a reference when determining the level of cut for the selected Calcar Body. The Osteotomy Guide has two notches that align with the level of the 45mm and 55mm Calcar Body build-up lengths. These notches are referenced off the tip of the greater trochanter. The guide also has markings to align the +0mm head center with each body build-up height. These markings are also referenced off the tip of the greater trochanter (Fig. 7).

The Cone and Calcar Bodies are designed for use when defects extend below the lesser trochanter or proximal bone stock is of poor quality. The Cone and Calcar Bodies are designed for use when defects extend below the lesser trochanter.

WARNING: The standard junction ZMR Hip System should only be used when full proximal support will be achieved in the area of the plasma spray. This is necessary because without full proximal support, the mid-stem junction is vulnerable to fracture. If such proximal support cannot be achieved, evaluate the use of the ZMR XL.

To prepare the proximal femur for the conical portion of any proximal body style, choose a Distal Pilot that matches the size of the distal implant that will be used. For example, if the distal femoral canal was prepared for a 15mm distal stem, select the size 15 Distal Pilot. The Distal Pilot diameter is 14.3mm, or 0.7mm smaller than the implant diameter to allow for clearance down the canal while maintaining a secure fit. If, however, the diaphysis is overreamed (e.g., for a bowed stem), choose a Distal Pilot that more closely matches the final reamed size as the pilot is used to help center the Porous Body Conical Reamer.

Thread the **proper** Distal Pilot onto a Porous Body Conical Reamer that is 1 or 2 sizes below the anticipated body size and tighten it with the Distal Pilot Wrench (Fig. 8). The instrument tray contains three of each size Distal Pilot so they can be attached in advance to incrementally sized Porous Body Conical Reamers. Sequentially ream to the desired size proximal body and appropriate depth.


Fig. 7


Fig. 8

Apply lateral and slight posterior pressure on the reamer to help maintain a centered orientation when reaming. The three most proximal marks on the Porous Body Conical Reamer match the 35mm, 45mm, and 55mm body build-up heights of the implants and are referenced off the tip of the greater trochanter (Fig. 9). If using a Calcar Body, advance the Porous Body Conical Reamer only to the most distal line on the reamer. This line should be referenced to the level of the osteotomy (Fig. 9a). These body lengths are referenced off the +0mm femoral head center.


Fig. 9

Fig. 9a

To help avoid disassociation of the Distal Pilot from the reamer, do not run the Porous Body Conical Reamer in reverse. (If the reamer is used in reverse and the Distal Pilot disassociates from the reamer, use the threaded end of the Distal Pilot Wrench to retrieve the Distal Pilot.) To facilitate removal of the Distal Pilot from the Porous Body Conical Reamer, attach the T-handle to the reamer and loosen with the Distal Pilot Wrench. If using the Cone or Calcar Body, preparation of the proximal femur is now complete.

If using the Spout Body, preparation of the medial calcar area must be completed. Select the appropriate Distal Pilot and thread it onto the appropriate size Spout Mill Guide (Fig. 10). Tighten it with the Distal Pilot Wrench. The final Porous Body Conical Reamer used will dictate which Spout Mill Guide to select and the final Distal Reamer size used will dictate the appropriate Distal Pilot (Table 2).

Table 2: Selection of Spout Mill Components

Conical Reamer Selected	Spout Mill Assembly	
	Spout Mill Guide Required	Spout Mill Cutter Required
AA/AA-	AA/AA-	AA/AA-
A	A	A-B
B	B	A-B
C	C	C-D
D	D	C-D
E	E	E-F
F	F	E-F

Fig. 10

Pass the Spout Mill Guide/Distal Pilot assembly down the femoral canal to the appropriate level using the markings on the lateral side of the guide. The lines on the Spout Mill Guide match the 35mm, 45mm, and 55mm body heights and align with the top of the greater trochanter (Fig. 11). Insert the Spout Mill Guide to the same depth as the last Porous Body Conical Reamer used. For example, if the reamer reached the 45mm depth mark, then the 45mm depth mark on the Spout Mill Guide should be aligned with the top of the greater trochanter. The assembly should seat without force. Use the Stem Impactor (00-9965-040-02) to control rotation of the Spout Mill Guide during insertion. This will determine the final anteversion of the femoral prosthesis. It is important to maintain the proper rotational orientation of the Spout Mill Guide in order to achieve proper femoral anteversion.

Attach the Spout Mill Cutter to a power drill/reamer. Align the two round projections on the mill cutter with the top slot of the Spout Mill Guide. Check the version of the Spout Mill Guide by moving the Spout Mill Cutter down the guide and assessing the position of the cutter relative to the proximal medial calcar (Fig. 12).


Fig. 11


Fig. 12

When the orientation of the Spout Mill Guide is satisfactory, turn on the drill/driver. Pass the Spout Mill Cutter down the femur while maintaining lateral pressure on the cutter. The slots in the Spout Mill Guide will direct the cutter in a path that matches the geometry of the proximal body spout (Fig. 13). When the cutter reaches the bottom of the guide, turn off the drill/driver power and remove the Spout Mill Assembly. Do not remove the cutter while the power is on as this may remove additional bone. Preparation of the proximal femur is now complete.

Note: If preparation was made for a 45mm or 55mm Spout Body, a Cone Body of the same size can still be used if it is determined that more version is needed. Full proximal support in the aes of plasma spray. The use of the Spout Mill Guide may not be necessary due to existing bone loss.


Fig. 13

3 Trial Reduction

Before inserting the final prosthesis, perform a trial reduction to assess leg length, abductor muscle tension, joint stability, and range of motion. If it is difficult to seat the provisional in the femur, there may be a need for additional reaming and/or milling to remove impinging bone. Assess the fit again with the provisional construct. The Stem Provisionals are undersized one millimeter compared to the actual implants of the same size, so it is not possible to judge the ease of insertion of the definitive implant with the provisionals.

There are Distal Stem Provisionals and Proximal Body Provisionals that match all the implants in the ZMR Hip System. Based on the instruments used to prepare the femur, choose the corresponding provisional body and stem.

If a straight distal stem component will be used, attach the Proximal Body Provisional to the straight Distal Stem Provisional using the captured Compression Nut. While securing the Proximal Body Provisional with the Proximal Body Wrench, use the Nut Driver and Torque Wrench to torque the nut clockwise to 15N-m (130 in.-lbs.) (Fig. 14). Attach the preoperatively planned Femoral Head Provisional to the trunnion of the Proximal Body Provisional. Insert the provisional construct into the femoral canal and perform a trial reduction. To evaluate joint stability, leg length, and offset, insert different provisional femoral head/neck combinations and body options, if needed.

With the provisional construct in place, assess the ability of the femur to support the proximal body above the mid-stem modular junction. If the proximal body cannot be adequately supported, an alternative implant may be appropriate.


Fig. 14

If a bowed distal stem component will be used, follow the same procedure as above, with one exception. In order to orient the proximal body in the correct version relative to the bow of the distal stem, the provisional body and stem should be loosely secured prior to insertion down the femoral canal. Once the provisional construct is placed in the femur, adjust the Proximal Body Provisional to the correct orientation in relation to the Distal Stem Provisional. Use the Stem Alignment Guide) to check the orientation of the bow by rotating the arm to center over the neck and locking it in position (Figs. 15, 15a, & 15b).


Fig. 15a


Fig. 15


Fig. 15b

Note: Arrow on Knurled Knob must coincide with stem tip direction. This is opposite from the apex of the bone.

When this position is acceptable, tightly secure the Proximal Body Provisional to the Distal Stem Provisional by using the Proximal Body Wrench, Nut Driver, and Torque Wrench. Tighten the Compression Nut to 15N-m (130 in.-lbs.) while the provisional assembly is in the femur (Fig. 16).

Using the Provisional Body Wrench will counter the tightening torque and minimize the risk of fracturing the femur during in situ tightening. Once tightened, proceed with the trial reduction.

Mark the orientation of the body on the bone as a reference to achieve the same version with the final implant. Then remove the provisional components and note the position of the proximal body relative to the bow of the distal stem.


Fig. 16

4 Implantation

A – Back Table Assembly Technique

If desired, the implant can be completely assembled before insertion. This back-table assembly is typically used for a straight stem, but can also be used for a bowed stem if the orientation of the distal stem and proximal body is certain.

Before assembling the implant on the back table, inspect the body and stem tapers to make sure they are clean and dry.

With a bowed stem, it is particularly important to orient the proximal body component in the correct position relative to the distal stem component. The short axis of the oval indent on top of the stem is aligned with the bow. Place the proximal body component loosely onto the distal stem component so the alignment matches that of the provisional construct used earlier. Insert the Stem Alignment Guide into the stem and rotate the stem until proper version is achieved.

Press the body and stem components together by hand for initial locking. Insert the Junction Assembly Instrument through the body counterbore, and secure it by engaging the threads of the distal stem and turning the knurled handle of the Junction Assembly Instrument.

If the Junction Assembly Instrument does not easily thread onto the stem, do not force it. Instead, remove the Junction Assembly Instrument and realign it to prevent cross threading.

Once the assembly instrument is securely fastened, confirm the rotation between the proximal body and distal stem. Squeeze the two handles together until the needle deflects to the indicator mark on the handle (Fig. 17). If the handles squeeze completely together before the needle is adequately deflected, tighten the knurled handle more to spread the handles further apart. Then reapply pressure to the handles until the needle deflects to the indicator mark.


Fig. 17

Note: Incorrect orientation of the knurled knob arrow will result in incorrect implant assembly.

After the mid-stem modular junction has been securely locked, insert and tighten the Compression Nut, which is packaged with the distal stem. Use the Proximal Body Wrench to secure the proximal body, and insert the Torque Wrench with Nut Driver into the body counterbore to engage the Compression Nut. Turn the Torque Wrench clockwise to 15N-m (130 in.-lbs.) (Fig. 18). Do not overtighten the Compression Nut as this could compromise its function. The assembled construct is now ready for insertion down the femoral canal.

Fig. 18

Insert the construct down the femoral canal by hand until it will no longer advance. Assess the implant for proper rotational alignment before impacting (Fig. 19). Then insert the Stem Impactor into the counterbore of the body, aligning it with the oval slot of the distal stem. Use the Mallet to drive the Stem Impactor and seat the implant in the femur.

The markings on the impactor serve as a reference (keyed off the greater trochanter) to the corresponding body height to assess the proper implant depth (Fig. 20). If the implant does not advance with each blow of the Mallet, stop insertion and remove the component. Then ream or mill additional bone from the areas that are preventing insertion and insert the component again.

With the implant construct in place, assess the support of the proximal body to ensure full proximal support has been achieved.

WARNING: The standard junction ZMR Hip System should only be used when full proximal support will be achieved in the area of the plasma spray. This is necessary because without full proximal support, the mid-stem junction is vulnerable to fracture. If such proximal support cannot be achieved, evaluate the use of the ZMR XL.

Fig. 19


Fig. 20

B – “Loose” Assembly Technique

If the alignment of the proximal body and the distal stem is uncertain, assemble the components loosely. Do not engage the taper, but loosely thread the Compression Nut, which is packaged with the distal stem, onto the distal stem.

Insert the loosely assembled prosthesis into the canal, allowing the proximal body to rotate into the appropriate orientation relative to the distal stem. Use the Stem Impactor and Mallet to impact the prosthesis to its final position. The markings on the impactor serve as a reference (keyed off the greater trochanter) to the corresponding body height to assess the proper implant depth (Fig. 20). The implant should advance with each moderate blow of the Mallet. If it does not, remove the implant and perform additional reaming.

Rotate the proximal body into the desired version. Then use the Nut Driver to remove the Compression Nut. Insert the Junction Assembly Instrument through the body counterbore and secure it by engaging the threads of the distal stem and turning the knurled handle.

When securely fastened, squeeze the two handles together until the needle deflects to the indicator mark on the handle. If this is not achieved on the first attempt, retighten the assembly and repeat this step.

Remove the Junction Assembly Instrument. Then re-insert and tighten the Compression Nut. Use the Proximal Body Wrench to secure the proximal body, and insert the Torque Wrench with Nut Driver into the body counterbore to engage the nut. Turn the Torque Wrench clockwise to 15N-m. Do not overtighten the Compression Nut as this could compromise its function.

With the implant construct in place, assess the support of the proximal body to ensure full proximal support has been achieved.

WARNING: The standard junction ZMR Hip System should only be used when full proximal support will be achieved in the area of the plasma spray. This is necessary because without full proximal support, the mid-stem junction is vulnerable to fracture. If such proximal support cannot be achieved, evaluate the use of the ZMR XL.


Fig. 21

5 Femoral Head Assembly

Once the implant is fully seated in the femoral canal, place the selected Femoral Head Provisional onto the taper of the implant. Perform a trial reduction to assess joint stability, range of motion, and restoration of leg length and offset.

When the appropriate femoral head implant is confirmed, remove the Femoral Head Provisional. Thoroughly clean and dry the taper on the femoral head component. Then place the selected femoral head on the taper and twist it slightly to secure it. **Impact it by striking the Head Impactor with one sharp blow of the Mallet.** Test the security of the head fixation by trying to remove it by hand.

Reduce the hip and assess leg length, range of motion, stability, and abductor tension.

Note: Always check that the neck taper and inside taper of the femoral head are clean and dry before impaction. Also, do not impact the femoral head onto the taper before driving the prosthesis down the femoral canal as the femoral head may loosen during impaction of the implant.

6 Wound Closure

After obtaining hemostasis, insert a Hemovac® Wound Drainage Device, if desired. Then close the wound in layers.

Postoperative Management

The postoperative management of patients with the ZMR Porous Revision Implant is determined by the surgical technique, patient bone quality, patient activity level, fit of the implant, and the surgeon's judgment. Weight bearing after revision surgery requires more external support for a longer period due to the nature of the extensive surgery and bone disruption.

Note: Accepted practices of postoperative care should be followed. The patient must be informed and made aware of the limitations of total joint reconstruction and the necessity of limiting weight and physical activity to protect the femoral stem from unnecessary stresses. See package insert for full patient counseling information.

Implant Removal

In the event that the ZMR Porous Revision Implant must be removed, specially designed instruments are available to remove the implant. Make sure that all bone and soft tissue proximal to the implant are removed. Begin by removing the Compression Nut from the implant. It is important to use the Proximal Body Wrench to apply a counter-rotating force to stabilize the implant and not stress the femur when using the Torque Wrench and Nut Driver to loosen the nut (Fig. 22). Failure to do so may result in fracture of the femur. Once the Proximal Body Wrench has secured the implant body, attach the Nut Driver to the Torque Wrench and insert it into the body counterbore to engage the Compression Nut. Then turn the Torque Wrench counterclockwise to loosen the nut.

Fig. 22

Remove the Compression Nut and attach the Distal Stem Extractor to the distal stem through the counterbore in the body by threading the extractor onto the threads in a clockwise direction. Engage as many threads as possible until tight (Fig. 23). **Be careful not to cross-thread the Distal Stem Extractor and stem.** If the Distal Stem Extractor is difficult to turn while threading it onto the stem, remove and realign the extractor.


Fig. 23

When the Distal Stem Extractor is completely engaged, attach the Slaphammer to the extractor by turning the Slaphammer handle clockwise to engage the threads of the extractor. Turn the handle until tight. Reverse impact the stem from the canal by moving the Slaphammer weight (Fig. 24).

The leg should be firmly held and the direction of force of the Slaphammer should be in line with the implant and femur. Extracting in a direction that is not in line with the femur may cause fracture of the femur.

Fig. 24

An additional option using an extractor hook is available for use only with the Spout or Calcar Bodies. The Spout Bodies and Calcar Bodies have extraction holes to accept the Extractor Hook (Fig. 25). Thread the Extractor Hook to the Slaphammer by turning the hook in a clockwise direction. Then reverse impact the stem from the canal by moving the Slaphammer weight.

Note: If the Compression Nut cannot be removed, impact a carbide punch into the inferior surface of the implant neck and drive the implant out of the canal with a Mallet.



Fig. 25

What ZMR Instruments Do I Need?

When not using Over-the-Junction Instrumentation:


COMPONENT	KIT DESCRIPTION	KIT NUMBER
Cone Body	Porous Revision General Instrument Set	00-9975-000-45
	Cone Body Provisional Instrument Set	00-9970-000-02
	Porous Proximal Preparation Instrument Set #1	00-9975-000-39
Calcar Body	Porous Revision General Instrument Set	00-9975-000-45
	Calcar Body Provisional Instrument Set	00-9970-000-03
	Porous Proximal Preparation Instrument Set #1	00-9975-000-39
Spout Body	Porous Revision General Instrument Set	00-9975-000-45
	Spout Body Provisional Instrument Set	00-9970-000-01
	Porous Proximal Preparation Instrument Set #2	00-9975-000-40
Straight Porous & Splined Stems	Porous Revision General Instrument Set	00-9975-000-45
	Straight Stem Provisional Instrument Set	00-9961-000-03
	Straight IM Reamer 10mm – 19.5mm Instrument Set	00-9975-000-06
	Straight IM Reamer 20mm – 26mm Instrument Set	00-9975-000-07
Bowed Porous, Porous Slotted & Splined Stems	Porous Revision General Instrument Set	00-9975-000-45
	Bowed Stem Provisional Instrument Set	00-9961-000-02
	Pressure Sentinel Intramedullary Reamers Instrument Set	00-9975-000-11
Tapered Distal Stem	Porous Revision General Instrument Set	00-9975-000-45
	Tapered “Crossover” Instrument Set	00-9965-000-02
	Taper Provisional Instrument Set	00-9972-000-00

1 General Instrument Kit


Prod. No	Description
00-9975-000-45	Porous Revision General Instrument Kit
00-9976-025-00	Case
00-9975-055-00	Tray 1
00-9976-026-00	Base
00-9975-099-00	Lid
00-7895-028-01	28MM Femoral Head Provisional (-3.5)
00-7895-028-02	28MM Femoral Head Provisional (0)
00-7895-028-03	28MM Femoral Head Provisional (+3.5)
00-7803-028-14	28MM MIS Slotted Femoral Head Provisional (+7)
00-7895-028-05	28MM Femoral Head Provisional (+10.5)
00-7895-032-01	32MM Femoral Head Provisional (-3.5)
00-7895-032-02	32MM Femoral Head Provisional (0)
00-7895-032-03	32MM Femoral Head Provisional (+3.5)
00-7803-032-14	32MM MIS Slotted Femoral Head Provisional (+7)
00-7895-032-05	32MM Femoral Head Provisional (+10.5)
00-9975-031-02	Distal Stem Extractor
00-6601-002-00	Extractor Hook
00-9027-058-00	Femoral Head Impactor
00-6551-006-00	Slaphammer
00-9975-035-02	Nut Driver
00-9975-033-02	Junction Assembly Pliers
00-9975-036-00	Socket
00-9975-036-10	Proximal Body Wrench
00-9965-040-02	Provisional Stem Impactor
00-9975-038-00	T- Handle
00-9975-041-00	Torque Wrench
00-9965-042-00	Osteotomy Guide
00-9955-071-03	Loose Assembly Sleeve
00-9965-140-00	Stem Alignment Guide

2 Spout Body Provisional Instrument Kit


Prod. No	Description
00-9970-000-01	Spout Body Provisional Instrument Kit
00-9976-040-00	Case
00-9976-042-00	Tray 1
00-9976-043-00	Tray 2
00-9976-044-00	Tray 3
00-9975-099-00	Lid
00-9930-016-36	AA-Body, Reduced Offset 35mm Build-Up (36 x 35 Neck)
00-9970-016-36	AA-Body, Reduced Offset 35mm Build-Up (36 x 35 Neck)
00-9930-016-40	AA-Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9970-016-40	AA-Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9970-017-40	A Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9970-018-40	B Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9970-019-40	C Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9970-020-40	D Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9970-021-40	E Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9970-023-40	F Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9977-017-45	A Body, Standard Offset 45mm Build-Up (40 x 45 Neck)
00-9977-018-45	B Body, Standard Offset 45mm Build-Up (40 x 45 Neck)
00-9977-019-45	C Body, Standard Offset 45mm Build-Up (40 x 45 Neck)
00-9977-017-55	A Body, Standard Offset 55mm Build-Up (40 x 55 Neck)
00-9977-018-55	B Body, Standard Offset 55mm Build-Up (40 x 55 Neck)
00-9977-019-55	C Body, Standard Offset 55mm Build-Up (40 x 55 Neck)
00-9970-017-46	A Body, Extended Offset 35mm Build-Up (46 x 35 Neck)
00-9970-018-46	B Body, Extended Offset 35mm Build-Up (46 x 35 Neck)
00-9970-019-46	C Body, Extended Offset 35mm Build-Up (46 x 35 Neck)
00-9970-020-46	D Body, Extended Offset 35mm Build-Up (46 x 35 Neck)
00-9970-021-46	E Body, Extended Offset 35mm Build-Up (46 x 35 Neck)
00-9970-023-46	F Body, Extended Offset 35mm Build-Up (46 x 35 Neck)
00-9971-017-45	A Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9971-018-45	B Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9971-019-45	C Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9971-020-45	D Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9971-021-45	E Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9971-023-45	F Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9971-017-55	A Body, Extended Offset 55mm Build-Up (46 x 55 Neck)
00-9971-018-55	B Body, Extended Offset 55mm Build-Up (46 x 55 Neck)
00-9971-019-55	C Body, Extended Offset 55mm Build-Up (46 x 55 Neck)
00-9971-020-55	D Body, Extended Offset 55mm Build-Up (46 x 55 Neck)
00-9971-021-55	E Body, Extended Offset 55mm Build-Up (46 x 55 Neck)
00-9971-023-55	F Body, Extended Offset 55mm Build-Up (46 x 55 Neck)

3 Cone Body Provisional Instrument Kit


Prod. No	Description
00-9970-000-02	Cone Body Provisional Instrument Kit
00-9976-035-00	Case
00-9976-037-00	Tray 1
00-9976-038-00	Tray 2
00-9976-039-00	Tray 3
00-5967-018-00	Lid
00-9932-016-35	AA-Body, Reduced Offset 35mm Build-Up (36 x 35 Neck)
00-9931-016-35	AA-Body, Reduced Offset 35mm Build-Up (36 x 35 Neck)
00-9933-016-35	AA-Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9976-016-35	AA-Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9976-017-35	A Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9976-018-35	B Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9976-019-35	C Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9976-020-35	D Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9976-021-35	E Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9976-023-35	F Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9973-017-35	A Body, Extended Offset 35mm Build-Up (46 x 35 Neck)
00-9973-018-35	B Body, Extended Offset 35mm Build-Up (46 x 35 Neck)
00-9973-019-35	C Body, Extended Offset 35mm Build-Up (46 x 35 Neck)
00-9973-020-35	D Body, Extended Offset 35mm Build-Up (46 x 35 Neck)
00-9973-021-35	E Body, Extended Offset 35mm Build-Up (46 x 35 Neck)
00-9973-023-35	F Body, Extended Offset 35mm Build-Up (46 x 35 Neck)
00-9976-017-45	A Body, Standard offset 45mm Build-Up (40 x 45 Neck)
00-9976-018-45	B Body, Standard offset 45mm Build-Up (40 x 45 Neck)
00-9976-019-45	C Body, Standard offset 45mm Build-Up (40 x 45 Neck)
00-9976-017-55	A Body, Standard Offset 55mm Build-Up (40 x 55 Neck)
00-9976-018-55	B Body, Standard Offset 55mm Build-Up (40 x 55 Neck)
00-9976-019-55	C Body, Standard Offset 55mm Build-Up (40 x 55 Neck)
00-9973-017-45	A Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9973-018-45	B Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9973-019-45	C Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9973-020-45	D Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9973-021-45	E Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9973-023-45	F Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9973-017-55	A Body, Extended Offset 55mm Build-Up (46 x 55 Neck)
00-9973-018-55	B Body, Extended Offset 55mm Build-Up (46 x 55 Neck)
00-9973-019-55	C Body, Extended Offset 55mm Build-Up (46 x 55 Neck)
00-9973-020-55	D Body, Extended Offset 55mm Build-Up (46 x 55 Neck)
00-9973-021-55	E Body, Extended Offset 55mm Build-Up (46 x 55 Neck)
00-9973-023-55	F Body, Extended Offset 55mm Build-Up (46 x 55 Neck)

4 Calcar Body Provisional Instrument Kit


Prod. No	Description
00-9970-000-03	Calcar Body Provisional Instrument Kit
00-9978-030-00	Case
00-9976-032-00	Tray 1
00-9976-033-00	Tray 2
00-9975-099-00	Lid
00-9937-016-35	AA-Body, Reduced Offset 35mm Build-Up (36 x 35 Neck)
00-9936-016-35	AA-Body, Reduced Offset 35mm Build-Up (36 x 35 Neck)
00-9938-016-35	AA-Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9978-016-35	AA-Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9978-017-35	A Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9978-018-35	B Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9978-019-35	C Body, Standard Offset 35mm Build-Up (40 x 35 Neck)
00-9978-017-45	A Body, Standard Offset 45mm Build-Up (40 x 45 Neck)
00-9978-018-45	B Body, Standard Offset 45mm Build-Up (40 x 45 Neck)
00-9978-019-45	C Body, Standard Offset 45mm Build-Up (40 x 45 Neck)
00-9978-017-55	A Body, Standard Offset 55mm Build-Up (40 x 55 Neck)
00-9978-018-55	B Body, Standard Offset 55mm Build-Up (40 x 55 Neck)
00-9978-019-55	C Body, Standard Offset 55mm Build-Up (40 x 55 Neck)
00-9974-017-45	A Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9974-018-45	B Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9974-019-45	C Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9974-020-45	D Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9974-021-45	E Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9974-023-45	F Body, Extended Offset 45mm Build-Up (46 x 45 Neck)
00-9974-017-55	A Body, Extended Offset 55mm Build-Up (46 x 55 Neck)
00-9974-018-55	B Body, Extended Offset 55mm Build-Up (46 x 55 Neck)
00-9974-019-55	C Body, Extended Offset 55mm Build-Up (46 x 55 Neck)
00-9974-020-55	D Body, Extended Offset 55mm Build-Up (46 x 55 Neck)
00-9974-021-55	E Body, Extended Offset 55mm Build-Up (46 x 55 Neck)
00-9974-023-55	F Body, Extended Offset 55mm Build-Up (46 x 55 Neck)

5 Straight Stem Provisional Instrument Kit



Prod. No	Description
00-9961-000-03	Straight Stem Provisional Instrument Kit
00-9976-045-00	Case
00-9976-046-00	Base
00-9975-099-00	Lid
00-9961-120-21	12.0mm x 115mm
00-9961-120-22	12.0mm x 170mm
00-9961-135-21	13.5mm x 115mm
00-9961-135-22	13.5mm x 170mm
00-9961-150-21	15.0mm x 115mm
00-9961-150-22	15.0mm x 170mm
00-9961-165-21	16.5mm x 115mm
00-9961-165-22	16.5mm x 170mm
00-9961-180-21	18.0mm x 115mm
00-9961-180-22	18.0mm x 170mm
00-9961-195-21	19.5mm x 115mm
00-9961-195-22	19.5mm x 170mm
00-9961-210-22	21.0mm x 170mm
00-9961-225-22	22.5mm x 170mm
00-9961-240-22	24.0mm x 170mm
00-9961-255-22	25.5mm x 170mm

6 Bowed Stem Provisional Instrument Kit



Prod. No	Description
00-9961-000-02	Bowed Stem Provisional Instrument Kit
00-9976-050-00	Case
00-9976-052-00	Tray 1
00-9976-053-00	Tray 2
00-9976-051-00	Base
00-9975-099-00	Lid
00-9961-120-32	12.0mm x 170mm
00-9961-135-32	13.5mm x 170mm
00-9961-150-32	15.0mm x 170mm
00-9961-165-32	16.5mm x 170mm
00-9961-180-32	18.0mm x 170mm
00-9961-195-32	19.5mm x 170mm
00-9961-210-32	21.0mm x 170mm
00-9961-225-32	22.5mm x 170mm
00-9961-240-32	24.0mm x 170mm
00-9961-255-32	25.5mm x 170mm
00-9961-135-33	13.5mm x 220mm
00-9961-150-33	15.0mm x 220mm
00-9961-165-33	16.5mm x 220mm
00-9961-180-33	18.0mm x 220mm
00-9961-195-33	19.5mm x 220mm
00-9961-210-33	21.0mm x 220mm
00-9961-225-33	22.5mm x 220mm
00-9961-240-33	24.0mm x 220mm
00-9961-255-33	25.5mm x 220mm
00-9963-150-34	15.0mm x 260mm
00-9963-165-34	16.5mm x 260mm
00-9963-180-34	18.0mm x 260mm
00-9963-195-34	19.5mm x 260mm
00-9963-210-34	21.0mm x 260mm
00-9963-225-34	22.5mm x 260mm
00-9963-240-34	24.0mm x 260mm
00-9963-255-34	25.5mm x 260mm

7 Porous Proximal Preparation Instrument Kit #1



Prod. No	Description
00-9975-000-39	Porous Proximal Preparation Instrument Kit
00-9976-055-00	Case
00-9976-057-00	Tray 1
00-9976-056-00	Base
00-9975-099-00	Lid
00-9975-001-16	AA/AA-Body, Conical Reamer
00-9975-001-17	A Body, Conical Reamer
00-9975-001-18	B Body, Conical Reamer
00-9975-001-19	C Body, Conical Reamer
00-9975-001-20	D Body, Conical Reamer
00-9975-001-21	E Body, Conical Reamer
00-9975-001-23	F Body, Conical Reamer
00-9975-002-12	12.0 Distal Pilots (2)
00-9975-002-13	13.5 Distal Pilots (2)
00-9975-002-15	15.0 Distal Pilots (2)
00-9975-002-16	16.5 Distal Pilots (2)
00-9975-002-18	18.0 Distal Pilots (2)
00-9975-002-19	19.5 Distal Pilots (2)
00-9975-002-21	21.0 Distal Pilots (2)
00-9975-002-22	22.5 Distal Pilots (2)
00-9975-002-24	24.0 Distal Pilots (2)
00-9975-002-25	25.5 Distal Pilots (2)
00-9975-039-00	Distal Pilot Wrench

8 Porous Proximal Preparation Instrument Kit #2



Prod. No	Description
00-9975-000-40	Porous Proximal Preparation Instrument Kit
00-9976-060-00	Case
00-9976-062-00	Tray 1
00-9976-061-00	Base
00-9975-099-00	Lid
00-9975-002-12	12.0 Distal Pilot
00-9975-002-13	13.5 Distal Pilot
00-9975-002-15	15.0 Distal Pilot
00-9975-002-16	16.5 Distal Pilot
00-9975-002-18	18.0 Distal Pilot
00-9975-002-19	19.5 Distal Pilot
00-9975-002-21	21.0 Distal Pilot
00-9975-002-22	22.5 Distal Pilot
00-9975-002-24	24.0 Distal Pilot
00-9975-002-25	25.5 Distal Pilot
00-9975-003-16	AA/AA-Body, Spout Mill Cutter
00-9975-003-17	A-B Body, Spout Mill Cutter
00-9975-003-19	C-D Body, Spout Mill Cutter
00-9975-003-21	E-F Body, Spout Mill Cutter
00-9975-005-16	AA/AA-Body, Spout Mill Guide
00-9975-005-17	A Body, Spout Mill Guide
00-9975-005-18	B Body, Spout Mill Guide
00-9975-005-19	C Body, Spout Mill Guide
00-9975-005-20	D Body, Spout Mill Guide
00-9975-005-21	E Body, Spout Mill Guide
00-9975-005-23	F Body, Spout Mill Guide

9 Straight IM Reamer 10 – 19.5mm Instrument Kit


Prod. No	Description
00-9975-000-06	Straight IM Reamer 10 – 19.5mm Instrument Kit
00-9976-065-00	Case
00-9976-067-00	Tray 1
00-9976-066-00	Base
00-9975-099-00	Lid
00-7891-010-00	10.0mm Straight IM Reamer
00-7891-010-05	10.50mm Straight IM Reamer
00-7891-011-00	11.00mm Straight IM Reamer
00-7891-011-05	11.50mm Straight IM Reamer
00-7891-012-00	12.00mm Straight IM Reamer
00-7891-012-05	12.50mm Straight IM Reamer
00-7891-013-00	13.00mm Straight IM Reamer
00-7891-013-05	13.50mm Straight IM Reamer
00-7891-014-00	14.00mm Straight IM Reamer
00-7891-014-05	14.50mm Straight IM Reamer
00-7891-015-00	15.00mm Straight IM Reamer
00-7891-015-05	15.5mm Straight IM Reamer
00-7891-016-00	16.0mm Straight IM Reamer
00-7891-016-05	16.5mm Straight IM Reamer
00-7891-017-00	17.0mm Straight IM Reamer
00-7891-017-05	17.5mm Straight IM Reamer
00-7891-018-00	18.0mm Straight IM Reamer
00-7891-018-05	18.5mm Straight IM Reamer
00-7891-019-00	19.0mm Straight IM Reamer
00-7891-019-05	19.5mm Straight IM Reamer

10 Straight IM Reamer 20 – 26mm Instrument Kit


Prod. No	Description
00-9975-000-07	Straight IM Reamer 20mm – 26mm Instrument Kit
00-9976-015-07	Case
00-9976-017-00	Tray 1
00-9976-016-00	Base
00-9975-099-00	Lid
00-7891-020-00	20.0mm Straight IM Reamer
00-7891-020-05	20.5mm Straight IM Reamer
00-7891-021-00	21.0 mm Straight IM Reamer
00-7891-021-05	21.5mm Straight IM Reamer
00-7891-022-00	22.0mm Straight IM Reamer
00-7891-022-05	22.5mm Straight IM Reamer
00-7891-023-00	23.0mm Straight IM Reamer
00-7891-023-05	23.5mm Straight IM Reamer
00-7891-024-00	24.0mm Straight IM Reamer
00-7891-024-05	24.5mm Straight IM Reamer
00-7891-025-00	25.0mm Straight IM Reamer
00-7891-025-05	25.5mm Straight IM Reamer
00-7891-026-00	26.0mm Straight IM Reamer

11 Taper Provisional Instrument Kit


Prod. No	Description
00-9972-000-00	Taper Provisional Instrument Kit
00-9976-016-00	Case
00-9975-591-00	Tray 1
00-9976-090-00	Base
00-9975-099-00	Lid
00-9975-040-02	Stem Impactor
00-9962-014-13	14mm x 135mm Revision Taper Stem Provisional
00-9962-015-13	15mm x 135mm Revision Taper Stem Provisional
00-9962-016-13	16mm x 135mm Revision Taper Stem Provisional
00-9962-017-13	17mm x 135mm Revision Taper Stem Provisional
00-9962-018-13	18mm x 135mm Revision Taper Stem Provisional
00-9962-019-13	19mm x 135mm Revision Taper Stem Provisional
00-9962-014-18	14mm x 185mm Revision Taper Stem Provisional
00-9962-015-18	15mm x 185mm Revision Taper Stem Provisional
00-9962-016-18	16mm x 185mm Revision Taper Stem Provisional
00-9962-017-18	17mm x 185mm Revision Taper Stem Provisional
00-9962-018-18	18mm x 185mm Revision Taper Stem Provisional
00-9962-019-18	19mm x 185mm Revision Taper Stem Provisional
00-9962-020-18	20mm x 185mm Revision Taper Stem Provisional
00-9962-022-18	22mm x 185mm Revision Taper Stem Provisional
00-9962-016-23	16mm x 235mm Revision Taper Stem Provisional
00-9962-017-23	17mm x 235mm Revision Taper Stem Provisional
00-9962-018-23	18mm x 235mm Revision Taper Stem Provisional
00-9962-019-23	19mm x 235mm Revision Taper Stem Provisional
00-9962-020-23	20mm x 235mm Revision Taper Stem Provisional
00-9962-022-23	22mm x 235mm Revision Taper Stem Provisional
00-7895-028-01	28mm Femoral Head Provisional (-3.5)
00-7895-028-02	28mm Femoral Head Provisional (0)
00-7895-028-03	28mm Femoral Head Provisional (+3.5)
00-7803-028-14	28mm MIS Slotted Femoral Head Provisional (+7)
00-7895-028-05	28mm Femoral Head Provisional (+10.5)
00-7895-032-01	32mm Femoral Head Provisional (-3.5)
00-7895-032-02	32mm Femoral Head Provisional (0)
00-7895-032-03	32mm Femoral Head Provisional (+3.5)
00-7803-032-14	32mm MIS Slotted Femoral Head Provisional (+7)
00-7895-032-05	32mm Femoral Head Provisional (+10.5)

12 Pressure Sentinel Intramedullary Reamers Instrument Kit


Prod. No	Description	Prod. No	Description
00-9975-000-11	Pressure Sentinel Intramedullary Reamers Instrument Kit	00-2228-017-05	17.5mm Flex Reamer
00-2228-040-00	Case	00-2228-018-00	18.0mm Flex Reamer
00-2228-042-00	Tray 1	00-2228-018-05	18.5mm Flex Reamer
00-2228-043-00	Tray 2	00-2228-019-00	19.0mm Flex Reamer
00-2228-044-00	Tray 3	00-2228-019-05	19.5mm Flex Reamer
00-2228-045-00	Tray 4	00-2228-020-00	20.0mm Flex Reamer
00-2228-041-00	Base	00-2228-020-05	20.5mm Flex Reamer
00-9975-099-00	Lid	00-2228-021-00	21.0mm Flex Reamer
00-2228-008-00	8.0mm Flex Reamer	00-2228-021-05	21.5mm Flex Reamer
00-2228-008-05	8.5mm Flex Reamer	00-2228-022-00	22.0mm Flex Reamer
00-2228-009-00	9.0mm Flex Reamer	00-2228-022-05	22.5mm Flex Reamer
00-2228-009-05	9.5mm Flex Reamer	00-2228-023-00	23.0mm Flex Reamer
00-2228-010-00	10.0mm Flex Reamer	00-2228-023-05	23.5mm Flex Reamer
00-2228-010-05	10.5mm Flex Reamer	00-2228-024-00	24.0mm Flex Reamer
00-2228-011-00	11.0mm Flex Reamer	00-2228-024-05	24.5mm Flex Reamer
00-2228-011-05	11.5mm Flex Reamer	00-2228-025-00	25.0mm Flex Reamer
00-2228-012-00	12.0mm Flex Reamer	00-2228-025-05	25.5mm Flex Reamer
00-2228-012-05	12.5mm Flex Reamer	00-2228-026-00	26.0mm Flex Reamer
00-2228-013-00	13.0mm Flex Reamer	00-2228-026-05	26.5mm Flex Reamer
00-2228-013-05	13.5mm Flex Reamer	00-2228-027-00	27.0mm Flex Reamer
00-2228-014-00	14.0mm Flex Reamer	00-2228-030-00	T Handle Extractor
00-2228-014-05	14.5mm Flex Reamer	00-5044-012-00	Universal Adapters (2)
00-2228-015-00	15.0mm Flex Reamer	00-9965-081-10	Diameter Gauge
00-2228-015-05	15.5mm Flex Reamer		
00-2228-016-00	16.0mm Flex Reamer		
00-2228-016-05	16.5mm Flex Reamer		
00-2228-017-00	17.0mm Flex Reamer		

Spout Body Options

Size (mm)	Size w/Plasma (mm)	Build-up Options (mm)	Offset Options (mm)	Body Length Options (mm)
AA-	16.7	35	36	75
	16.7	35	40	75
AA	16.7	35	36	80
	16.7	35	40	80
A	17.7	35*	40*, 46	80
	17.7	45	46	90
	17.7	55	46	100
B	18.9	35*	40*, 46	80
	18.9	45	46	90
	18.9	55	46	100
C	20	35*	40*, 46	80
	20	45	46	90
	20	55	46	100
D	21.2	35*	40*, 46	80
	21.2	45	46	90
	21.2	55	46	100
E	22.5	35*	40*, 46	80
	22.5	45	46	90
	22.5	55	46	100
F	23.7	35*	40*, 46	80
	23.7	45	46	90
	23.7	55	46	100


* 40mm offset only provided with 35mm build-up Spout Body. Build-up, offset, and body length dimensions shown at +0mm head center.

Cone Body Options

Size (mm)	Size w/Plasma (mm)	Build-up Options (mm)	Offset Options (mm)	Body Length Options (mm)
AA-	16.7	35	36	75
	16.7	35	40	75
AA	16.7	35	36	80
	16.7	35	40	80
A	17.7	35	40, 46	80
	17.7	45	40, 46	90
	17.7	55	40, 46	100
B	18.9	35	40, 46	80
	18.9	45	40, 46	90
	18.9	55	40, 46	100
C	20	35	40, 46	80
	20	45	40, 46	90
	20	55	40, 46	100
D	21.2	35	40, 46	80
	21.2	45	40, 46	90
	21.2	55	40, 46	100
E	22.5	35	40, 46	80
	22.5	45	40, 46	90
	22.5	55	40, 46	100
F	23.7	35	40, 46	80
	23.7	45	40, 46	90
	23.7	55	40, 46	100


Build-up, offset, and body length dimensions shown at +0mm head center.

Calcar Body Options

Size (mm)	Size w/Plasma (mm)	Build-up Options (mm)	Offset Options (mm)	Body Length Options (mm)
AA-	16.7	35	36	75
	16.7	35	40	75
AA	16.7	35	36	80
	16.7	35	40	80
A	17.7	35	40, 46	80
	17.7	45	40, 46	90
	17.7	55	40, 46	100
B	18.9	35	40, 46	80
	18.9	45	40, 46	90
	18.9	55	40, 46	100
C	20	35	40, 46	80
	20	45	40, 46	90
	20	55	40, 46	100
D	21.2	45	46	90
	21.2	55	46	100
E	22.5	45	46	90
	22.5	55	46	100
F	23.7	45	46	90
	23.7	55	46	100


Build-up, offset, and body length dimensions shown at +0mm head center.

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S denotes available in a straight stem, B denotes available in a bowed stem, L denotes available in slotted bowed stem.

Assembled Implant Length*

Body Length (mm)	+	Stem Length (mm)	=	Implant Length (mm)
75		115		190
75		170		245
75		220		295
75		160		335
80		115		195
80		170		250
80		220		300
80		260		340
90		115		205
90		170		260
90		220		310
90		260		350
100		115		215
100		170		270
100		220		320
100		260		360

* Using +0mm Head Center

References

1. Agnew SG, Stewart CL, Skinner R, et al. Design influences on pressure and temperature generation from intramedullary reamers in vitro and in vivo models. Presented at the Orthopaedic Trauma Association Annual Meeting, Tampa, FL; September, 1995.

This documentation is intended exclusively for physicians and is not intended for laypersons. Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendations. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice in whole or in part. Please refer to the package inserts for important product information, including, but not limited to, contraindications, warnings, precautions, and adverse effects.

Contact your Zimmer representative or visit us at www.zimmer.com

The CE mark is valid only if it is also printed on the product label.

+H12497999010200/\$110831111-

97-9990-102-00 1108-H05 7.5ML Printed in USA ©2011 Zimmer, Inc.

ZMR®
**Over-the-Junction
Instruments for
Revision Hip
Arthroplasty**
Surgical Technique

Introduction

The *ZMR* Over-the-Junction (OTJ) Instrumentation includes four key components; 1) OTJ Proximal Reamers 2) OTJ Taper Protectors 3) OTJ Trochanteric Starter Reamers and 4) Metal Proximal Provisionals. These instruments can be used with any combination of *ZMR* standard junction or *ZMR* XL junction implants. *ZMR* Over-the-Junction Instrumentation is designed to increase the predictability of otherwise unpredictable femoral revision hip surgery.

OTJ Proximal Reamers

The OTJ Proximal Reamers are designed to provide a consistent and reproducible mechanism for preparing the proximal femur and may be used in conjunction with either the distal provisionals or the final distal stem implant. They are matched for each proximal body size for both *ZMR* standard and *ZMR* XL junction stems. Visualization windows in the cutting flute area of each reamer allow the surgeon to accurately determine complete reaming depth.

OTJ Taper Protectors

ZMR OTJ Instrumentation includes taper protectors that are unique to the *ZMR* OTJ system. Not only do these unique instruments completely protect the locking taper junction from inadvertent reamer strikes or damage from bone fragments, they also serve as a solid and safe alignment guide for the proximal reamers to follow.

OTJ Trochanteric Starter Reamers

It is not uncommon for the greater trochanter to protrude over the femoral canal, which could impede proper placement of the distal implant/Taper Protector assembly or the distal provisional and Taper Protector. For this reason, Trochanteric Starter reamers are an integral part of the OTJ Instrumentation, providing the surgeon with an effective tool to use when removing trochanteric or metaphyseal bone.

Metal Proximal Provisionals

ZMR metal proximal provisionals are radiographically visible and are designed to safely lock onto the definitive distal implant. They resist twisting on the distal implant, and can be easily removed or repositioned if necessary. These trials help the surgeon determine the proper implant version, build-up, and offset to effectively maximize joint stability and to optimize the biomechanical restoration of the hip joint function.

WARNING: The standard junction *ZMR* Hip System should only be used when full proximal support will be achieved in the area of the plasma spray. This is necessary because without full proximal support, the mid-stem junction is vulnerable to fracture. If such proximal support cannot be achieved, evaluate the use of the *ZMR* XL.

1 Distal Preparation

The ZMR OTJ reamers can be used after the distal trial stem or definitive distal implant has been fully seated in the femur. As such, follow the appropriate ZMR surgical technique from the beginning through the desired point when proximal preparation is initiated.


Refer to the current surgical techniques as indicated –

- For Porous, Porous Slotted and Splined Stems: *ZMR Porous Revision Hip Prosthesis Surgical Technique (97-9990-002-00)*
- For Tapered Distal Stems: *ZMR Crossover Instruments and Surgical Technique (97-9990-007-00)*
- For XL Junction Implants: *ZMR XL Taper and Porous Hip Prostheses Surgical Technique (97-9921-002-00)*.

When impacting the distal trial/stem, be sure to reference the depth marks on the stem inserter (# 00-9955-040-02) to the tip of the greater trochanter as these correspond to the +0mm head center of the various proximal body heights and will aid in selecting the appropriate proximal body.

Note: If selecting the ZMR Spout Bodies, OTJ reaming must be initiated using the distal trial as milling must occur prior to implantation of the definitive stem implant.

WARNING: The standard junction ZMR Hip System should only be used when full proximal support will be achieved in the area of the plasma spray. This is necessary because without full proximal support, the mid-stem junction is vulnerable to fracture. If such proximal support cannot be achieved, evaluate the use of the ZMR XL.

2 Using the Trochanteric Starter Reamer

Ensure there is adequate space to allow the Taper Protector to be attached to the distal stem or trial. Use the Trochanteric Starter Reamer to clear away any part of the greater trochanter and metaphyseal bone that impedes this process. Be sure to use only the standard size Trochanteric Starter Reamer (Item # 00-9975-007-04) for the standard junction stems/trials, and the XL size (Item # 00-9975-007-05) for the XL stems/trials. If preparing for a tapered distal stem, leave the last distal reamer in place and position the Trochanteric Starter Reamer over the distal reamer's shaft. This will function as a guide. When using a porous distal stem, attach the reamer to power, start rotation of the reamer, then gently move the reamer down the axis of the femoral canal to remove the impeding bone. Stop advancing the reamer when the markings on the shaft align with the tip of the greater trochanter.

Note: If using the Trochanteric Starter Reamer after the insertion of the distal implant, care must be taken to avoid reamer contact with the distal stem taper as this could compromise the integrity of the locking junction.

3 Assembling the Taper Protector

Select the appropriately sized Taper Protector (Standard Junction Item #00-9975-007-00 / XL Junction Item #00-9975-007-01) and thread onto the proximal threads of the distal stem implant or the distal provisional. If unable to do so by hand, use the Stem Inserter (Item #00-9955-040-02) by placing it into the proximal inserter feature on the Taper Protector and turn in a clockwise motion until it stops.

Note: Do not over-tighten. If the Taper Protector will not fully thread onto the distal stem/provisional or stops advancing, remove the Taper Protector. Visually inspect all of the threads of the protector and stem / provisional to ensure there is no debris or damage and re-attach.

Note: Do not use the Taper Protector to impact or further seat the distal stem as damage to the threads or the distal stem may result.

4 Proximal Reaming

Select a Proximal Reamer that matches the taper junction size (Standard or XL) and Taper Protector. Surgeons have the option of either starting one size smaller than the templated size of the proximal body implant to be used (Example- The templated body size is “C”; start reaming with the size “B” reamer) or proceeding directly to the templated size. This is surgeon preference and may be based upon familiarization with the instruments, surgical technique, and their assessment of the patient’s bone quality. Position the reamer over the Taper Protector and on power with gentle pressure, insert the reamer until it fully seats to the hard-stop provided by the Taper Protector. To prevent the risk of the reamer binding on the greater trochanter, be sure the reamer is rotating before engaging the femur. Windows in the reamer flutes allow the surgeon to visualize when the reamer has fully bottomed out against the Taper Protector.

Warning: Only ream over a properly attached Taper Protector. Reaming without the Taper Protector may cause damage to the implant’s locking taper and will compromise implant function and durability. Do not use the XL OTJ Proximal Reamer in combination with the standard size Taper Protector. This mismatch in sizing could compromise prepared fit between the proximal femur and the proximal body as well as potentially damage the instruments.

If choosing to implant a Cone Body, the preparation of the proximal femur is now complete. If selecting to use a Calcar or Spout body, please follow the preparation steps outlined in the ZMR Porous Revision Hip Prosthesis Surgical Technique.

Taper Protector is visible

5 Remove Taper Protector

Remove the Taper Protector from the distal implant either by hand or by inserting the Stem Inserter (Item #00-9955-040-02) into the proximal inserter feature and turning counter clockwise.

6 Using Proximal Provisionals for Determining Version, Offset, and Build-Up

The surgeon has two main proximal body provisional options, and selection depends on whether a distal provisional was inserted and used for proximal preparation or the actual distal implant was placed and used.

- If using the distal provisional, the non-metallic proximal provisionals must be used. For the Cone Body, select the appropriate body size from the Cone Body Provisional Instrument Kit #00-9970-000-02. For the Calcar Body, use Calcar Body Provisional Instrument Kit #00-9970-000-03 and for the Spout Bodies, use Spout Body Provisional Instrument Kit #00-9970-000-01. In all instances, select the proximal provisional which corresponds to the body size of the last proximal reamer used. Follow the steps outlined in the appropriate surgical technique and finish the surgery.
- If the definitive distal implant was used during proximal preparation, the surgeon may either use the non-metallic proximal provisionals as indicated above, or use the Metal Proximal Provisionals found in the ZMR OTJ Instrumentation Kit #KT-9975-000-00. Important Note: The Metal Proximal Provisionals only replicate the head center height (build-up) and offset. The bodies are smaller in diameter than the actual implants, and therefore should not be used to assess the proximal fill of the final implant. If the metal provisionals are to be used on the final distal stem, proceed to Step #7.

Proximal body
metal provisional


7

Using Metal Proximal Provisionals

Select the desired Metal Proximal Provisional (MPP), identify the desired proximal version, and place the provisional on the distal implant. Place the Stem Inserter in the inserter feature located on the proximal lateral aspect of the provisional neck and, using only the heel of the hand, tap onto the stem. Do not strike with a mallet. Perform trial reduction.

Note: The MPP locks onto the taper of the definitive distal stem. As such, there may be a small mark left on the proximal stem taper. This will NOT affect performance of the stem taper as the proximal implant body locks in an area distal to this mark.

Mark the desired stem version on the femur with a bovie or marking pen. Remove the MPP by inserting the Stem Inserter into the screw located within the proximal shoulder of the provisional and turn counter clockwise one-quarter to one-half turn. If the provisional becomes locked on, the provisional does have an extraction hole that can be used with the Extractor hook (Item #00-6601-002-00) which is found in the Porous Revision General Instrument Kit #00-9975-000-45.

Used to remove MPP from stem by turning counter clockwise

Used to secure MPP on stem


8

Implant Assembly

Assemble the proximal body implant and finish the procedure per the instructions in the appropriate surgical technique.

Ordering Information

Instrument Scope List

ZMR Over-the-Junction Instrument Kit

Order: KT-9975-000-00 (Includes 1 each of the following)

Item Number	Description	Item Number	Description	Item Number	Description
1	00-9975-007-16 Proximal Reamer – AA Body	10	00-9975-007-01 Taper Protector - XL	20	00-9996-046-45 Proximal Provisional 46 x 45 STD
2	00-9975-007-17 Proximal Reamer – A Body	11	00-9975-007-04 Trochanteric Starter Reamer - STD	21	00-9996-046-55 Proximal Provisional 46 x 55 STD
3	00-9975-007-18 Proximal Reamer – B Body	12	00-9975-007-05 Trochanteric Starter Reamer - XL	22	00-9924-040-78 Proximal Provisional 40 x 78 XL
4	00-9975-007-19 Proximal Reamer – C Body	13	00-9996-036-30 Proximal Provisional 36 x AA- STD	23	00-9924-040-83 Proximal Provisional 40 x 83 XL
5	00-9975-007-20 Proximal Reamer – D Body	14	00-9996-036-35 Proximal Provisional 36 x 35 STD	24	00-9924-040-93 Proximal Provisional 40 x 93 XL
6	00-9975-007-21 Proximal Reamer – E Body	15	00-9996-040-30 Proximal Provisional 40 x AA- STD	25	00-9924-040-13 Proximal Provisional 40 x 103 XL
7	00-9975-007-23 Proximal Reamer – F Body	16	00-9996-040-35 Proximal Provisional 40 x 35 STD	26	00-9924-045-78 Proximal Provisional 45 x 78 XL
8	00-9975-007-02 Proximal Reamer – XL Body	17	00-9996-040-45 Proximal Provisional 40 x 45 STD	27	00-9924-045-83 Proximal Provisional 45 x 83 XL
9	00-9975-007-00 Taper Protector - STD	18	00-9996-040-55 Proximal Provisional 40 x 55 STD	28	00-9924-045-93 Proximal Provisional 45 x 93 XL
		19	00-9996-046-35 Proximal Provisional 46 x 35 STD	29	00-9924-045-13 Proximal Provisional 45 x 103 XL

This documentation is intended exclusively for physicians and is not intended for laypersons. Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendations. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice in whole or in part. Please refer to the package inserts for important product information, including, but not limited to, contraindications, warnings, precautions, and adverse effects.

Contact your Zimmer representative or visit us at www.zimmer.com

The CE mark is valid only if it is also printed on the product label.

+H12497999010900/\$110831111-

ZMR® Crossover Instruments

Abbreviated
Surgical Technique

Introduction

ZMR Crossover Instruments facilitate the combination of any Porous Proximal Body with any Taper Stem in the ZMR System, thus allowing extensive fixation in the femur. The instrumentation and implant combinations provide the opportunity to treat wide variances in patient anatomy.

Tapered Stems

ZMR Revision Taper Components were designed to achieve secure distal fixation in the femur using a sharply splined and tapered distal stem with a roughened titanium surface. The tapered distal stem is designed to wedge into the femoral medullary canal, transferring axial and bending forces, while the splines are press-fit into the bone to provide rotational stability. The roughened *Tivanium*® Ti-6Al-4V Alloy surface allows bone ongrowth.^{1,2,3} A bevel at the distal end of the stem is a design feature intended to increase the ease of insertion, to better accommodate the bow of the femur, and decrease the potential for distal femoral cortical perforation.

Published clinical results of other stems using this design philosophy (tapered, splined with a roughened titanium alloy surface) in femoral revision surgery have been impressive. These results show the favorable remodeling of proximal femoral bone stock when excessive bone loss was present.^{4,5}


Porous Body Options

Multiple sizes in each body type allow for metaphyseal filling, proximal fixation, and proximal support of the prosthesis.

Spout Body (A) – The Spout Body helps achieve medial fill in the femur, contributing to initial rotational stability.

Cone Body (B) – The Cone Body addresses cases where the Spout Body's medial geometry is not desired. The Cone Body provides the opportunity for infinite version adjustment.

Calcar Body (C) – The Calcar Body has a medial collar to help resist subsidence.

Simplicity

Crossover Instruments, consisting of proximal and distal reamers, are contained in just one additional tray, facilitating ease of the surgical procedure. Cannulated proximal reamers match the porous body geometry and facilitate a similar technique compared to the existing Revision Taper System.

Please note:


1. The Distal Crossover Reamer is identifiable by its titanium nitrided (gold color) cutting flutes.
2. The Proximal Crossover Reamer is discernable from the standard Porous Proximal Reamer by its titanium nitrided (gold color) shaft.

Distal Canal Preparation


1 Ream Distal Femur

Based on preoperative templating, ream the distal femur to the appropriate size and depth. Hand reaming is recommended.

2 Final Ream


Advance the final reamer until it corresponds to one of the three body heights (35mm, 45mm, or 55mm) and leave the reamer in place. Remove the T-handle from the reamer and note the stem length chosen. Consider obtaining a cross-table A/P radiograph to confirm proper sizing and positioning in the femur.

Proximal Femoral Preparation

3 Assess Need for Sleeve


From the stem length selected, determine if a sleeve is required and if needed, place the appropriate sleeve(s) on the distal reamer.

4 Ream Proximal Femur


Starting with a cannulated proximal cone reamer one or two sizes smaller than the templated size, ream the proximal femur to the appropriate size. Match the depth of the reamer to the depth noted in Step One.

WARNING: The standard junction ZMR Hip System should only be used when full proximal support will be achieved in the area of the plasma spray. This is necessary because without full proximal support, the mid-stem junction is vulnerable to fracture. If such proximal support cannot be achieved, evaluate the use of the ZMR XL.


5 Prepare Femur for Spout Body (Optional)

If a Spout Body is to be used in the patient, the **Porous Body Conical Reamer** (non-cannulated) must be used to ream the shelf for the Spout Mill Guide.

5a Mill Femur (Optional)

The Spout Mill Guide and Cutter can then be used to prepare the medial metaphysis. The 12mm Distal Pilot should be used to center the Guide in the canal.

Trial Reduction

6 Assemble Trial

Assemble the appropriate proximal body and distal stem provisionals and position the distal tip bevel anteriorly.

7 Check Positioning in Femur

Insert the assembled trial into the proximal femur after visualizing the medullary canal to retrieve any debris from reaming. If trial will not completely seat, check position of anterior bevel and/or re-ream distal canal, if necessary.

8 Lock Assembled Trial

Once assembled, tighten the Compression Nut with the Torque Wrench Assembly. Perform trial reduction.

10 Lock Assembled Construct

Insert the Compression Nut and tighten to 15N-m.


Note: Do not overtighten the Compression Nut as this could compromise its function.

Implant Insertion


9 Assemble Implant

Assemble the final implants with the Junction Assembly Tool while taking care to replicate the rotational orientation of the anterior bevel.


11 Insert Implant

Insert the assembled implant to the appropriate depth while maintaining proper anteversion of the femoral neck.

Ordering Information

Assembled Implant Length*

Body Length (mm)	+	Stem Length (mm)	=	Implant Length (mm)
75		135		210
75		185		260
75		235		310
80		135		215
80		185		265
80		235		315
90		135		225
90		185		275
90		235		325
100		135		235
100		185		285
100		235		335

* Using +0mm Head Center

Note: For more details, refer to the ZMR Hip System brochure (97-9990-001-00), or the Revision Taper Surgical Technique (97-9982-002-00) and Porous Revision Surgical Technique (97-9990-002-00).

References

1. Feighan JE, Goldberg VM, Davy D, et al. The influence of surface-blasting in the incorporation of titanium-alloy implants in a rabbit intramedullary model. *J Bone Joint Surg.* 1995;77-A;9:1380-1395.
2. Zweymuller KA, Lintner FK, Semlitsch MF. Biologic fixation of a press-fit titanium hip joint endoprosthesis. *Clin Orthop.* 1988;235:195-206.
3. Lintner F, Zweymuller KA, Brand G. Tissue reactions to titanium endo-prosthesis. *J Arthroplasty.* 1986;1;3:183-195.
4. Michelinakis E, Papapolychroniou T, Vafiadis J. The use of a cementless femoral component for the management of bone loss in revision hip arthroplasty. *Hosp for Joint Diseases.* 1996;55;1:28-32.
5. Hartwig CH, Bohm P, Czech U, et al. The Wagner revision stem in alloarthroplasty of the hip. *Arch Orthop Trauma Surg.* 1996;115:5-9.

Instruments

Prod. No.	Description
00-9965-000-02	Crossover Instrument Set
*00-9965-001-16	Proximal Crossover Reamer Size AA Body
00-9965-001-17	Proximal Crossover Reamer Size A Body
00-9965-001-18	Proximal Crossover Reamer Size B Body
00-9965-001-19	Proximal Crossover Reamer Size C Body
00-9965-001-20	Proximal Crossover Reamer Size D Body
00-9965-001-21	Proximal Crossover Reamer Size E Body
00-9965-001-23	Proximal Crossover Reamer Size F Body
00-9965-013-13	Distal Crossover Reamer Size 13
00-9965-014-13	Distal Crossover Reamer Size 14
00-9965-015-13	Distal Crossover Reamer Size 15
00-9965-016-13	Distal Crossover Reamer Size 16
00-9965-017-13	Distal Crossover Reamer Size 17
00-9965-018-13	Distal Crossover Reamer Size 18
00-9965-019-13	Distal Crossover Reamer Size 19
00-9965-020-18	Distal Crossover Reamer Size 20
00-9965-022-18	Distal Crossover Reamer Size 22
00-9975-037-05	Conical Reamer Sleeve
00-9976-010-00	Crossover Reamer Case

* Contact Zimmer Customer Service for availability.

This documentation is intended exclusively for physicians and is not intended for laypersons. Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendations. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice in whole or in part. Please refer to the package inserts for important product information, including, but not limited to, contraindications, warnings, precautions, and adverse effects.

Contact your Zimmer representative or visit us at www.zimmer.com

The CE mark is valid only if it is also printed on the product label.

87-6203-831-22

Rev. – September 2011

Printed in U.S.A.

© 2011, Zimmer, Inc.

ZMR[®] HIP SYSTEM

zimmer

Zimmer, Inc.
1800 West Center Street
Warsaw, Indiana 46580
USA
1-800-348-2759 (US only)
+1-574-372-4999

R_x only

Zimmer U.K. Ltd.
9 Lancaster Place
South Marston Park
Swindon, SN3 4FP, UK

87-6203-831-22

Rev. – September 2011
Printed in U.S.A.
© 2011, Zimmer, Inc.

ZMR® HIP SYSTEM

zimmer

Zimmer, Inc.
1800 West Center Street
Warsaw, Indiana 46580
USA
1-800-348-2759 (US only)
+1-574-372-4999

Rx only

Zimmer U.K. Ltd.
9 Lancaster Place
South Marston Park
Swindon, SN3 4FP, UK

ENGLISH

The CE mark is valid only if it is also printed on the product label.

ZMR® HIP SYSTEM

- Before using a product placed on the market by Zimmer, the operating surgeon should study carefully the following recommendations, warnings and instructions, as well as the available product-specific information (e.g., product literature, written surgical technique). Zimmer is not liable for complications that may arise from use of the device in circumstances outside of Zimmer's control including, but not limited to, product selection and deviations from the device's indicated uses or surgical technique.
- Due to the acquisition of pre-existing product lines, Zimmer has initiated a testing program to evaluate the compatibility of these devices with implants and components made or distributed by all Zimmer companies, which include Zimmer GmbH (previously Centerpulse Orthopedics Ltd.), Zimmer Trabecular Metal Technology, Inc. (previously Implex Corp.), Zimmer U.K. Ltd., and Zimmer Austin, Inc. (previously Centerpulse Orthopedics, Inc.). Only authorized combinations should be used. To determine whether these devices have been authorized for use in a proposed combination, please contact your Zimmer sales representative or visit the Zimmer website: www.productcompatibility.zimmer.com. A printout of the

website information can also be obtained by calling Zimmer, Inc. Customer Service, 1-800-348-2759 (U.S.) or the local international access code +1-574-372-4999 (outside the U.S.). Former Centerpulse and Implex products that are now packaged in Zimmer boxes, and for which compatibility could be an issue, have been labeled “former Centerpulse” and “former Implex” to provide clarification for the user.

DESCRIPTION

- The *ZMR* Hip System is a modular hip system for cementless revision hip arthroplasty. The stem assembly is manufactured from *Titanium*® Ti-6Al-4V Alloy and utilizes modular junctions between the head and neck, and the metaphyseal body and diaphyseal stem segments. A variety of body and stem components are provided to achieve fixation and restore joint kinematics.
- The head/neck modular connection of the femoral stem assembly is a 12/14 taper designed to mate with the corresponding femoral head component. The *ZMR* Hip Prosthesis may be used with wrought *Zimaloy*® Cobalt-Chromium-Molybdenum Alloy femoral heads or, where available and approved for use, the following ceramic heads:
 - *Zimmer*® Alumina Ceramic Femoral Heads (6418 and 6428 series)
 - *BIOLOX*† forte Ceramic Heads (12.28.05/07; 12.32.05/07)
 - *BIOLOX*® delta Ceramic Heads (8775 series)
 - *BIOLOX* OPTION Ceramic Head/Adapter System (8777 series)
 - *Cerasul*® Ceramic-on-Ceramic Tribological Solution Ceramic Heads (407.0101/03)
 - *Sulox*™ Ceramics Material Ceramic Heads (17.28.05/07; 17.32.05/07)

INDICATIONS

The *ZMR* Hip System is indicated for cementless revision hip arthroplasty. Arthroplasty should be performed only when more conservative methods of treatment have failed to provide symptomatic relief or when there is progressive disability.

INFORMATION FOR USE

Special instruments are available and necessary to assemble, disassemble, and accurately install the implant components.

CONTRAINDICATIONS

- Contraindications of the *ZMR* Hip System include:
 - Use of this device when a less invasive procedure would be sufficient
 - Skeletal immaturity
 - Loss of abductor musculature in the affected limb
 - Poor bone quantity or quality, where there is a lack of adequate bone to support the implant(s)
 - Poor skin coverage around the hip joint
 - Neuromuscular disease (e.g., Charcot’s joint) in the affected limb

- Additional contraindications for the *ZMR* Revision Taper and Revision Porous hip prostheses include:
 - Any clinical situation in which full proximal support cannot be achieved

WARNINGS

- The *ZMR* Revision Taper and Revision Porous hip prostheses should only be used when full proximal support will be achieved in the area of the plasma spray. This is necessary because without full proximal support, the mid-stem junction is vulnerable to fracture. If such proximal support cannot be achieved, evaluate the use of the *ZMR* XL hip prostheses.
- Proper function and longevity of the hip can be compromised if the patient engages in high levels of physical activity or does not comply with postoperative restrictions. Other medical disabilities that create or result in an unusual gait may exacerbate stresses on the femoral component.
- The load-bearing capacity of the implant can be compromised during surgery by notching, scratching, or striking the prosthesis, repeated assembly/disassembly of modular components, or failing to provide metaphyseal support to the implant.
- The *ZMR* Hip System is used for indications which generally require extensive surgery. This tends to render a prosthetic hip joint less stable than hip joints that follow conventional total hip replacement. There is, therefore, greater likelihood of postoperative dislocation. Postoperative positioning and range of motion activities must be carefully supervised.
- The smaller sized femoral components are designed for patients with a small intramedullary canal and/or metaphyseal region of the femur. The stem cross-sectional geometry has been decreased to accommodate this special anatomical requirement. This reduction in size results in a corresponding reduction in the fatigue strength of the implant. Because of this, the prospective patient should be screened carefully. These smaller femoral stems are suitable only for a patient who has a low-to-moderate activity level. Patients who are young, heavy and those who engage in high physical activity levels are not suitable candidates for the smaller sized femoral prostheses.
- Do not use any component if damage is found or caused during setup or insertion.
- Improper selection, placement, positioning, and fixation of the implant components may result in unusual stress conditions reducing the service life of the prosthetic implants.
- Avoid varus stem placement. Varus stem placement results in higher stresses on the stem, increasing the risk of fatigue fracture of the stem.
- Do not use forceful tapping to drive the femoral stem into place.
- Intraoperative radiographs to ascertain whether proximal femoral fracture has occurred may be prudent if:
 - Preparation of the medullary canal or insertion of the prosthesis is more difficult than predicted.
 - There is sudden loss of resistance when the broach or implant is impacted.
- This device is for single patient use only. Do not reuse.
- Do not use this product for other than labeled indications (off-label use).

PRECAUTIONS

- Continued surveillance for new or recurrent sources of infection should be continued as long as the device is in place.
- Do not mate the titanium alloy components of this system with stainless steel. Galvanic corrosion may develop with this combination.
- Do not assemble the mating components without ensuring that the surfaces are clean of blood or debris. Failure to examine and ensure no particular matter is present could result in inadequate seating of one component upon the other and subsequent disassembly of the mated components or fracture of the implant.
- Repeated assembly and disassembly of the junction (head/neck and body to stem) could compromise the critical locking action of the Morse-type tapers. Use the provisional components during trial reductions. Change the head/neck component only when clinically necessary.

ADVERSE EFFECTS

- The following adverse effects have been reported:
 - Peripheral neuropathies
 - Deep wound infections
 - Perforation of the acetabulum or femur
 - Heterotopic bone formation
 - Metal sensitivity
 - Inflammatory reactions and osteolysis
 - Fatigue fracture
 - Wear
 - Vascular complications
 - Trochanteric problems
 - Subclinical nerve damage
 - Corrosion of metal implants
 - Early or late loosening of components
 - Pelvic, femoral, or acetabular fractures
 - Disassembly of modular components
 - Dislocation and subluxation

STERILITY

These devices are provided sterile by gamma irradiation, indicated by the **STERILE** symbol on the labeling. They remain sterile as long as the package integrity has not been violated. Inspect each package prior to use and do not use the component if any seal or cavity is damaged or breached or if the expiration date has been exceeded. Once opened, the component must be used, discarded, or resterilized.

STERILIZATION INSTRUCTIONS

- These sterilization instructions are consistent with ANSI/AAMI/ISO standards and guidelines. They should be used for sterile items that were opened but unused.
- Do not reuse instruments or devices labeled for single use only.
- In the event of inadvertent loss of sterility while preparing for surgery, sterile packaged, all metal implants may be resterilized only once for immediate use. This is subject to the exceptions listed below.

- **DO NOT RESTERILIZE:**
 - Single use only components that have been contaminated with body fluids or debris or previously implanted.
 - Components with a packaging expiration date that has been exceeded.
- Do not use the original plastic cavities or lids for resterilization. Single devices should be packaged in a medical grade sterilization pouch or wrap which conforms to the recommended specifications for steam sterilization provided in the table below. Ensure that the pouch or wrap is large enough to contain the devices without stressing the seals or tearing.
- Rinse porous components to remove lint or debris (using USP purified water).
- Aggressive cleaning with detergents and brushes may damage special features of the implant, such as fiber metal pads or bead coatings. Also, certain detergents may be difficult to rinse off polymer items, especially those made of silicone rubber.
- Items made from titanium and titanium alloys can form oxide layers from steam boiler treatment chemicals or detergent residues. While these oxides are biocompatible, they can obliterate etchings and stampings.
- Modular implant components must be sterilized separately to minimize potential bioburden buildup in the dead space and expansion/contraction stresses.

Recommended Sterilization/Resterilization Specifications

Follow the sterilizer manufacturer's instructions for loading patterns and selection of sterilization parameters. Drying times vary according to load size and should be increased for larger loads.

Single, All Metal Implants

Steam Sterilization

Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
Gravity Displacement	121°C (250°F)	30 minutes	15 minutes
Gravity Displacement	132°C (270°F)	15 minutes	
UK Pre-vacuum/Pulsating Vacuum ¹	134°C (273°F)	3 minutes	
Pre-vacuum/Pulsating Vacuum	132°C (270°F)	4 minutes	

¹ This cycle is not for use in the United States.

Please contact *Zimmer* at the following number if you have additional questions. In the USA, call 1-800-348-2759. For calls outside the USA, call the local international access code +1-574-267-6131.

PATIENT COUNSELING INFORMATION

Complications and/or failure of prosthetic implants are more likely to occur in patients with unrealistic functional expectations, heavy patients, physically active patients, and/or with patients who fail to follow through with the required rehabilitation program. Physical activity or trauma can result in loosening, wear, and/or fracture of the implant. The patient must be counseled about the capabilities of the implant and the impact it will have on his or her lifestyle. The patient must be instructed about all postoperative restrictions, particularly those related to occupational and sports activities and about the possibility that the implant or its components may wear out, fail or need to be replaced. The implant may not last the rest of the patient's life, or any particular length of time. Because prosthetic implants are not as strong, reliable, or durable as natural, healthy tissues/bones, all such devices may need to be replaced at some point.

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†BIOLOX® is a trademark of CeramTec AG
†BILOX® , CeramTec AG'nin tescilli markasıdır