

ZAKTUALIZOWANA PILNA INFORMACJA DOTYCZĄCA BEZPIECZEŃSTWA

Wycofanie wstępnie załadowanych soczewek EyeCee ONE i EyeCee ONE Crystal

13 lipca 2023 r.

Szanowny Kliencie!

Firma NIDEK dobrowolnie wycofuje określone partie wstępnie załadowanych soczewek wewnątrzgałkowych EyeCee ONE i EyeCee ONE Crystal w krajach, do których dostarczono wadliwe soczewki wewnątrzgałkowe.

Załączamy pełną listę soczewek wewnątrzgałkowych, których to dotyczy (patrz Załącznik A).

Wycofanie jest wynikiem otrzymania przez firmy NIDEK oraz Bausch + Lomb zgłoszeń dotyczących podwyższonego ciśnienia śródgałkowego u pacjentów, którym wszczepiono jedną z tych soczewek wewnątrzgałkowych. W okresie kwarantanny firma NIDEK przeprowadziła wewnętrzne dochodzenie i ustaliła, że przyczyną jest konkretna partia materiału zastosowanego w procesie powlekania.

Na dyszy wstrzykiwacza zastosowano środek powlekający (PVP) w celu umożliwienia swobodnego wstrzyknięcia soczewki wewnątrzgałkowej. W wyniku dochodzenia firma NIDEK określiła KONKRETNE PARTIE wstrzykiwacza z materiałem PVP, w przypadku których dochodziło do fizycznego blokowania drogi odpływu cieczy wodnistej. Prawdopodobnie skutkuje to występowaniem nieprawidłowo podwyższonego ciśnienia śródgałkowego.

Zastosowanie tych produktów może stwarzać ryzyko utrzymującego się podwyższonego ciśnienia śródgałkowego, a pacjenci, którym wszczepiono te soczewki, powinni być ściśle monitorowani pod kątem poziomu ciśnienia śródgałkowego. Wszelkie przypadki wysokiego ciśnienia śródgałkowego wymagają natychmiastowej uwagi i leczenia chirurgicznego lub leczenia przez odpowiedniego pracownika opieki zdrowotnej.

Zdrowie i bezpieczeństwo każdego, kto korzysta z naszych produktów, jest naszym najwyższym priorytetem i traktujemy przypadki takie jak ten bardzo poważnie, dlatego podejmujemy to dobrowolne działanie.

Poinformowaliśmy o tym działaniu wszystkie właściwe organy regulacyjne. Niniejszy dokument zawiera szczegóły dotyczące tego wycofania, w związku z tym dalsze informacje nie będą rozpowszechniane.

Działania, które ma podjąć użytkownik

Zgodnie z dokumentacją przechowywaną w firmie NIDEK Państwa placówka otrzymała produkty objęte wycofaniem. Uprzejmie prosimy o podjęcie następujących kroków w odpowiedzi na niniejszą informację dotyczącą bezpieczeństwa:

- Przejrzeć zapasy wszystkich produktów, które są objęte wycofaniem.
- Odseparować w swojej placówce wszystkie produkty objęte wycofaniem.
(Szczegółowe informacje na temat produktów objętych wycofaniem znajdują się w „Załączniku A”)
- Autoryzowany dystrybutor firmy NIDEK skontaktuje się z Państwem w sprawie wypełnienia „Formularza weryfikacji klienta” i zwrotu produktów objętych wycofaniem. Do czasu dalszego kontaktu ze strony dystrybutora prosimy o odseparowanie i zachowanie wszystkich produktów objętych wycofaniem.
- Prosimy o przekazanie tej wiadomości wszystkim, którym mogli Państwo przekazać soczewki wewnętrzne objęte wycofaniem.
- Niniejszą pilną informację dotyczącą bezpieczeństwa należy przekazać wszystkim użytkownikom produktów, których dotyczy problem, oraz innym osobom w Państwa organizacji. Prosimy o zachowanie tych informacji przynajmniej do czasu zakończenia związanych z nimi czynności.

Informacje kontaktowe

Dziękujemy za cierpliwość i przepraszamy za wszelkie niedogodności związane z tą sprawą.

W przypadku pytań dotyczących tego problemu należy skontaktować się z lokalnym przedstawicielem handlowym:

FRANCJA

Dział chirurgiczny: service-clients.france@bausch.com

☎ +33 (0) 4 67 13 47 49

WŁOCHY

Dział chirurgiczny: servizio_clienti@bausch.com

☎ 02 91483851

HISZPANIA

Dział chirurgiczny: Servicio.Clientes@bausch.com

☎ +34 902 381 010

PORTUGALIA

Dział chirurgiczny: Servicio.Clientes@bausch.com

☎ +351 808 203 178

HOLANDIA, BELGIA

Dział chirurgiczny: csbenelux@bausch.com

☎ Holandia: +31 (0) 20 203 54 01; Belgia: +32 3 280 82 40.

IRLANDIA

Dział chirurgiczny: uksurgorders@bausch.com

☎ +44(0)208 781 0000

SZWECJA, NORWEGIA, DANIA, FINLANDIA, ISLANDIA

Dział chirurgiczny: customerservice.nordic@bausch.com

☎ Szwecja i Islandia: +46 08 616 9570; Norwegia: +47 800 104 40; Dania: +45 808 809 90;

Finlandia: +358 0800 118 011

POLSKA

Dział chirurgiczny: zamowienia.chirurgia@bausch.com

☎ +48 801 080 023

NIEMCY, AUSTRIA

Dział chirurgiczny: kundenservice@bausch.com

☎ Niemcy: +49 (0) 800 589 3114; Austria: +43 (0) 800 241015

Europa Środkowo-Wschodnia (Bułgaria, Cypr, Czechy, Estonia, Grecja, Chorwacja, Węgry, Litwa, Rumunia)

Dział chirurgiczny: ELC_EASTERN_EUROPE@bausch.com

☎ +31 206582728

Autoryzowany przedstawiciel NIDEK w krajach UE:

[Audrey guillou@nidek.fr](mailto:Audrey_guillou@nidek.fr)

☎ +33 (0) 1 49 80 97 97

Z poważaniem



Imię i nazwisko: Katsuaki Tohyama



Tytuł: starszy menedżer

Dział zapewnienia jakości

NIDEK

Zaktualizowana pilna informacja dotycząca bezpieczeństwa

Przykłady etykiet

Numer seryjny na opakowaniu zewnętrznym	Numer seryjny na opakowaniu sterylizowanym
 <p>The label on the external packaging is white with black and blue text. It features the model number 'SZ-1 +23. 5D' at the top. Below it is a blue-bordered box containing the serial number 'SN UJX05300' and a date '2027-11'. A barcode is positioned below the box. At the bottom, a long alphanumeric string '(01) 04987669585351 (17) 271100 (21) UJX05300' is printed.</p>	 <p>The label on the sterilized packaging is white with black and blue text. It includes the manufacturer's name 'NIDEK CO., LTD.' and address '34-14, Maehama, Hiroishi-cho, Gamagori, Aichi 443-0038, Japan'. A barcode is present, followed by the CE mark and '0344'. The model number 'SZ-1 BAUSCH + LOMB +23. 5D' is printed in large letters. Below it, dimensions are listed: '(φB) : 6. 0mm' and '(φT) : 13. 0mm'. A blue-bordered box contains the serial number 'SN UJX05300'. The 'EyeCee ONE preloaded' logo and the code '78400-P376-E0' are also visible.</p>

Na etykiecie opakowania zewnętrznego i opakowania sterylizowanego w celu identyfikacji produktu umieszczono numery seryjne. 2-cyfrowy prefiks oznacza rok i miesiąc produkcji produktu. Szczegółowe informacje znajdują się w tabeli na następnej stronie.

PILNA INFORMACJA DOTYCZĄCA BEZPIECZEŃSTWA

Rok produkcji	Miesiąc produkcji	Pierwsza cyfra	Druga cyfra	Rok produkcji	Miesiąc produkcji	Pierwsza cyfra	Druga cyfra
2019	5	P	F	2021	9	R	3
2019	6	P	G	2021	10	R	4
2019	7	P	H	2021	11	R	5
2019	8	P	J	2021	12	R	6
2019	9	P	K	2022	1	R	7
2019	10	P	L	2022	2	R	8
2019	11	P	M	2022	3	R	9
2019	12	P	N	2022	4	R	0
2020	1	R	A	2022	5	U	A
2020	2	R	B	2022	6	U	B
2020	3	R	C	2022	7	U	C
2020	4	R	E	2022	8	U	E
2020	5	R	F	2022	9	U	F
2020	6	R	G	2022	10	U	G
2020	7	R	H	2022	11	U	H
2020	8	R	J	2022	12	U	J
2020	9	R	K	2023	1	U	K
2020	10	R	L	2023	2	U	L
2020	11	R	M				
2020	12	R	N				
2021	1	R	P				
2021	2	R	R				
2021	3	R	T				
2021	4	R	U				
2021	5	R	W				
2021	6	R	Y				
2021	7	R	1				
2021	8	R	2				

WAŻNE: Partie produkcyjne z okresu od września 2021 r. do listopada 2022 r. podlegają wycofaniu (zaznaczone na żółto). Szczegółowe informacje na temat numerów seryjnych podlegających wycofaniu można znaleźć w „Załączniku A”.

(Przykład: Ponieważ produkt UJX05300 pochodzi z grudnia 2022 r., nie podlega wycofaniu z rynku) Przynajmniej problemem są produkty pochodzące z konkretnego okresu produkcji. Dlatego można używać produktów, które nie są objęte wycofaniem.

Zaktualizowany formularz potwierdzenia otrzymania pilnej informacji dotyczącej bezpieczeństwa

Niniejszym potwierdzamy otrzymanie wyżej wymienionego powiadomienia o wycofaniu z dnia 13 lipca 2023 r.

Prosimy o wypełnienie formularza i odesłanie zgodnie z instrukcją dystrybutora firmy NIDEK.

Po zwróceniu produktów podlegających wycofaniu otrzymają Państwo zwrot pieniędzy.

Uprzejmie informujemy, że w przypadkach, gdy zwrot produktu nie będzie zgodny z instrukcjami, na przykład zwrot bez tego formularza, płatność może nastąpić z opóźnieniem.

Szczegóły produktu:

Wstępnie załadowane soczewki wewnątrzgałkowe EyeCee ONE i EyeCee ONE crystal

Prosimy o zapoznanie się z poniższym oświadczeniem i potwierdzenie (X):

- Zapoznaliśmy się z załączoną informacją dotyczącą bezpieczeństwa i potwierdzamy otrzymanie ostrzeżenia.

Prosimy o zapoznanie się z poniższymi stwierdzeniami i potwierdzenie (X) tego, które dotyczy Państwa placówki:

- W naszych zapasach nie ma żadnych soczewek wewnątrzgałkowych, których dotyczy problem.
- Mamy w naszych zapasach soczewki wewnątrzgałkowe, których dotyczy problem. W przypadku zaznaczenia prosimy wypełnić poniższą tabelę.

Prosimy podać całkowitą liczbę soczewek wewnątrzgałkowych objętych wycofaniem zidentyfikowanych w Państwa placówce:

Nazwa produktu	Nr partii	Ilość zużytych	Ilość zwracanych

Aby uzyskać numer autoryzacji zwrotu towaru (RMA) i umówić się na odbiór zidentyfikowanego produktu, prosimy o kontakt z autoryzowanym dystrybutorem firmy NIDEK w Państwa kraju/regionie pod wskazanymi numerami kontaktowymi.

Niniejszym zaświadczam, że wymieniony produkt został odseparowany, aby uniemożliwić jego użycie, i oczekuje na odbiór.

Data

Imię i nazwisko (drukowanymi literami)

Numer konta

Podpis

Nazwa placówki

Numer telefonu

Prosimy o wypełnienie, podpisanie i odesłanie niniejszego formularza do autoryzowanego dystrybutora firmy NIDEK w Państwa kraju/regionie na adres e-mail wskazany w danych kontaktowych.

Umowa o zachowaniu poufności: Informacje zawarte w tej wiadomości są poufne i przeznaczone do użytku przez adresata podanego powyżej. Jeśli nie jesteś zamierzonym odbiorcą ani pracownikiem lub agentem odpowiedzialnym za dostarczenie tej wiadomości do zamierzonego odbiorcy, niniejszym informujemy, że jakiegokolwiek ujawnianie, kopiowanie, dystrybucja informacji lub podejmowanie jakichkolwiek działań w oparciu o nie jest surowo zabronione. W przypadku omyłkowego otrzymania tego pisma prosimy o natychmiastowe powiadomienie nas telefonicznie, aby umówić się na zwrot oryginalnego dokumentu.

Attachment A

Model	Lot No.
EyeCee One Preloaded +23.5	R3X33119
EyeCee One Preloaded +22.0	R4X03191
EyeCee One Preloaded +22.0	R4X35350
EyeCee One Preloaded +22.0	R4X35355
EyeCee One Preloaded +22.0	R4X35358
EyeCee One Preloaded +20.5	R3X22070
EyeCee One Preloaded +20.5	R3X22071
EyeCee One Preloaded +20.5	R3X22073
EyeCee One Preloaded +20.5	R3X22079
EyeCee One Preloaded +21.0	R3X24368
EyeCee One Preloaded +26.5	R3X32923
EyeCee One Preloaded +26.5	R3X32924
EyeCee One Preloaded +21.5	R3X33321
EyeCee One Preloaded +21.0	R3X46758
EyeCee One Preloaded +21.0	R3X46759
EyeCee One Preloaded +24.0	R3X41530
EyeCee One Preloaded +24.0	R3X41531
EyeCee One Preloaded +24.0	R3X41534
EyeCee One Preloaded +26.0	R4X23114
EyeCee One Preloaded +26.0	R4X23117
EyeCee One Preloaded +26.0	R4X23118
EyeCee One Preloaded +23.0	R4X25746
EyeCee One Preloaded +23.0	R4X25748
EyeCee One Preloaded +21.5	R4X28334
EyeCee One Preloaded +21.5	R4X40268
EyeCee One Preloaded +21.5	R4X40269
EyeCee One Preloaded +19.5	R4X44461
EyeCee One Preloaded +19.5	R5X23160
EyeCee One Preloaded +19.5	R5X23161
EyeCee One Preloaded +23.5	R5X26370
EyeCee One Preloaded +23.5	R5X26377
EyeCee One Preloaded +20.0	R3X25510
EyeCee One Preloaded +20.0	R3X25516
EyeCee One Preloaded +18.5	R3X42598
EyeCee One Preloaded +18.5	R3X42599
EyeCee One Preloaded +19.0	R3X54751
EyeCee One Preloaded +19.0	R3X54753
EyeCee One Preloaded +17.5	R3X57100

EyeCee One Preloaded +17.5	R3X57104
EyeCee One Preloaded +23.0	R4X04731
EyeCee One Preloaded +23.0	R4X04732
EyeCee One Preloaded +15.5	R4X28512
EyeCee One Preloaded +21.5	R4X30196
EyeCee One Preloaded +24.0	R4X35103
EyeCee One Preloaded +24.0	R4X35105
EyeCee One Preloaded +24.0	R4X35109
EyeCee One Preloaded +21.0	R4X37393
EyeCee One Preloaded +21.0	R4X37396
EyeCee One Preloaded +19.5	R5X04071
EyeCee One Preloaded +23.5	R5X07844
EyeCee One Preloaded +23.5	R5X07847
EyeCee One Preloaded +26.0	R5X07890
EyeCee One Preloaded +26.0	R5X07895
EyeCee One Preloaded +26.0	R5X07898
EyeCee One Preloaded +22.0	R5X09981
EyeCee One Preloaded +22.0	R5X09982
EyeCee One Preloaded +22.0	R5X09983
EyeCee One Preloaded +22.0	R5X09988
EyeCee One Preloaded +22.0	R5X09989
EyeCee One Preloaded +24.5	R5X20584
EyeCee One Preloaded +20.5	R5X23323
EyeCee One Preloaded +20.5	R5X23328
EyeCee One Preloaded +20.5	R5X23329
EyeCee One Preloaded +22.5	R3X16205
EyeCee One Preloaded +22.5	R3X16206
EyeCee One Preloaded +22.5	R3X16207
EyeCee One Preloaded +23.5	R3X33444
EyeCee One Preloaded +23.5	R3X33445
EyeCee One Preloaded +23.5	R3X33449
EyeCee One Preloaded +22.5	R3X42681
EyeCee One Preloaded +23.5	R3X47576
EyeCee One Preloaded +19.0	R4X07450
EyeCee One Preloaded +19.0	R4X07451
EyeCee One Preloaded +19.0	R4X07452
EyeCee One Preloaded +19.0	R4X07453
EyeCee One Preloaded +19.0	R4X07454
EyeCee One Preloaded +21.0	R4X12130

EyeCee One Preloaded +21.0	R4X12133
EyeCee One Preloaded +21.0	R4X12134
EyeCee One Preloaded +21.0	R4X12135
EyeCee One Preloaded +21.0	R4X12136
EyeCee One Preloaded +21.0	R4X12138
EyeCee One Preloaded +23.0	R4X22943
EyeCee One Preloaded +23.0	R4X22944
EyeCee One Preloaded +17.5	R4X27896
EyeCee One Preloaded +17.5	R4X27897
EyeCee One Preloaded +17.5	R4X27898
EyeCee One Preloaded +15.5	R4X28506
EyeCee One Preloaded +21.5	R4X30390
EyeCee One Preloaded +21.5	R4X30395
EyeCee One Preloaded +21.5	R4X30396
EyeCee One Preloaded +21.5	R4X30397
EyeCee One Preloaded +21.5	R4X30398
EyeCee One Preloaded +21.5	R4X30399
EyeCee One Preloaded +18.0	R4X32114
EyeCee One Preloaded +18.0	R4X32115
EyeCee One Preloaded +18.0	R4X32116
EyeCee One Preloaded +18.0	R4X32117
EyeCee One Preloaded +23.0	R4X35500
EyeCee One Preloaded +23.0	R4X35501
EyeCee One Preloaded +27.0	R4X45051
EyeCee One Preloaded +27.0	R4X45052
EyeCee One Preloaded +29.0	R4X47912
EyeCee One Preloaded +19.5	R5X04052
EyeCee One Preloaded +19.5	R5X04053
EyeCee One Preloaded +26.0	R5X04422
EyeCee One Preloaded +26.0	R5X04435
EyeCee One Preloaded +26.0	R5X04436
EyeCee One Preloaded +26.0	R5X04437
EyeCee One Preloaded +19.5	R5X04690
EyeCee One Preloaded +19.5	R5X04695
EyeCee One Preloaded +19.5	R5X04696
EyeCee One Preloaded +26.5	R5X17274
EyeCee One Preloaded +26.5	R5X17277
EyeCee One Preloaded +22.0	R5X21130
EyeCee One Preloaded +22.0	R5X21131

EyeCee One Preloaded +22.0	R5X21132
EyeCee One Preloaded +22.0	R5X21137
EyeCee One Preloaded +22.0	R5X21139
EyeCee One Preloaded +24.5	R5X21201
EyeCee One Preloaded +24.5	R5X21202
EyeCee One Preloaded +24.5	R5X21203
EyeCee One Preloaded +24.5	R5X21204
EyeCee One Preloaded +24.5	R5X21206
EyeCee One Preloaded +20.5	R5X23340
EyeCee One Preloaded +20.5	R5X23343
EyeCee One Preloaded +20.5	R5X23344
EyeCee One Preloaded +20.5	R5X23345
EyeCee One Preloaded +20.5	R5X26103
EyeCee One Preloaded +20.5	R5X26104
EyeCee One Preloaded +21.0	R3X55464
EyeCee One Preloaded +21.0	R4X17670
EyeCee One Preloaded +21.0	R4X17673
EyeCee One Preloaded +21.0	R4X17676
EyeCee One Preloaded +21.0	R4X17678
EyeCee One Preloaded +21.0	R4X17679
EyeCee One Preloaded +26.0	R5X07870
EyeCee One Preloaded +26.0	R5X07871
EyeCee One Preloaded +26.0	R5X07872
EyeCee One Preloaded +26.0	R5X07873
EyeCee One Preloaded +26.0	R5X07874
EyeCee One Preloaded +26.0	R5X09070
EyeCee One Preloaded +19.0	R4X47783
EyeCee One Preloaded +23.5	R3X21991
EyeCee One Preloaded +23.5	R3X21992
EyeCee One Preloaded +23.5	R3X21993
EyeCee One Preloaded +23.5	R3X21994
EyeCee One Preloaded +23.5	R3X21995
EyeCee One Preloaded +16.0	R3X39840
EyeCee One Preloaded +16.0	R3X39843
EyeCee One Preloaded +19.5	R3X42431
EyeCee One Preloaded +19.5	R3X42432
EyeCee One Preloaded +19.5	R3X42433
EyeCee One Preloaded +19.5	R3X42437
EyeCee One Preloaded +22.0	R3X45260

EyeCee One Preloaded +22.0	R3X45266
EyeCee One Preloaded +22.0	R3X45267
EyeCee One Preloaded +22.0	R3X45268
EyeCee One Preloaded +22.0	R3X45269
EyeCee One Preloaded +17.0	R3X47532
EyeCee One Preloaded +17.0	R3X47533
EyeCee One Preloaded +17.0	R3X47534
EyeCee One Preloaded +17.0	R3X47535
EyeCee One Preloaded +23.0	R3X55150
EyeCee One Preloaded +23.0	R3X55151
EyeCee One Preloaded +23.0	R3X55152
EyeCee One Preloaded +23.0	R3X55159
EyeCee One Preloaded +17.5	R3X57111
EyeCee One Preloaded +17.5	R3X57114
EyeCee One Preloaded +17.5	R3X57115
EyeCee One Preloaded +17.5	R3X57118
EyeCee One Preloaded +20.5	R4X10761
EyeCee One Preloaded +20.5	R4X10763
EyeCee One Preloaded +22.5	R4X12060
EyeCee One Preloaded +21.0	R4X18093
EyeCee One Preloaded +21.0	R4X18095
EyeCee One Preloaded +21.0	R4X18096
EyeCee One Preloaded +21.0	R4X18097
EyeCee One Preloaded +21.0	R4X18098
EyeCee One Preloaded +21.0	R4X18099
EyeCee One Preloaded +22.5	R4X20510
EyeCee One Preloaded +22.5	R4X20517
EyeCee One Preloaded +22.5	R4X20518
EyeCee One Preloaded +22.5	R4X20519
EyeCee One Preloaded +20.5	R4X23680
EyeCee One Preloaded +20.5	R4X23686
EyeCee One Preloaded +20.5	R4X23687
EyeCee One Preloaded +20.5	R4X23688
EyeCee One Preloaded +24.0	R4X23801
EyeCee One Preloaded +24.0	R4X23802
EyeCee One Preloaded +24.0	R4X23803
EyeCee One Preloaded +24.0	R4X23804
EyeCee One Preloaded +24.0	R4X23806
EyeCee One Preloaded +23.0	R4X34562

EyeCee One Preloaded +21.5	R4X37301
EyeCee One Preloaded +21.5	R4X37304
EyeCee One Preloaded +21.5	R4X37305
EyeCee One Preloaded +21.5	R4X37306
EyeCee One Preloaded +21.5	R4X37307
EyeCee One Preloaded +21.5	R4X37309
EyeCee One Preloaded +20.0	R4X37691
EyeCee One Preloaded +20.0	R4X37692
EyeCee One Preloaded +20.0	R4X37693
EyeCee One Preloaded +20.0	R4X37694
EyeCee One Preloaded +20.0	R4X37695
EyeCee One Preloaded +20.0	R4X37699
EyeCee One Preloaded +26.5	R5X04542
EyeCee One Preloaded +18.0	R5X06985
EyeCee One Preloaded +18.0	R5X06986
EyeCee One Preloaded +18.0	R5X06989
EyeCee One Preloaded +27.0	R5X15252
EyeCee One Preloaded +27.0	R5X15253
EyeCee One Preloaded +26.5	R5X17283
EyeCee One Preloaded +26.0	R5X18397
EyeCee One Preloaded +26.0	R5X18398
EyeCee One Preloaded +18.0	R5X18416
EyeCee One Preloaded +26.0	R6X08599
EyeCee One Preloaded +29.0	R3X33582
EyeCee One Preloaded +11.0	R4X20314
EyeCee One Preloaded +21.5	R4X28315
EyeCee One Preloaded +21.5	R3X41994
EyeCee One Preloaded +20.0	R4X44873
EyeCee One Preloaded +26.0	R6X39482
EyeCee One Preloaded +18.5	R6X08420
EyeCee One Preloaded +18.5	R6X08421
EyeCee One Preloaded +18.5	R6X08422
EyeCee One Preloaded +18.5	R6X08423
EyeCee One Preloaded +18.5	R6X08424
EyeCee One Preloaded +18.5	R6X08425
EyeCee One Preloaded +18.5	R6X08426
EyeCee One Preloaded +18.5	R6X08427
EyeCee One Preloaded +18.5	R6X08429
EyeCee One Preloaded +18.5	R6X41239

EyeCee One Preloaded +15.0	R5X01252
EyeCee One Preloaded +18.0	R5X18435
EyeCee One Preloaded +13.0	R5X49104
EyeCee One Preloaded +14.0	R5X49130
EyeCee One Preloaded +21.0	R6X02452
EyeCee One Preloaded +24.5	R6X03149
EyeCee One Preloaded +19.0	R6X08306
EyeCee One Preloaded +20.5	R6X08515
EyeCee One Preloaded +22.5	R6X09289
EyeCee One Preloaded +12.0	R6X14977
EyeCee One Preloaded +17.0	R6X15221
EyeCee One Preloaded +25.0	R6X18104
EyeCee One Preloaded +26.0	R6X19681
EyeCee One Preloaded +23.5	R6X21689
EyeCee One Preloaded +22.0	R6X22179
EyeCee One Preloaded +24.0	R6X24315
EyeCee One Preloaded +20.0	R6X26634
EyeCee One Preloaded +21.5	R6X29533
EyeCee One Preloaded +23.0	R6X31894
EyeCee One Preloaded +27.0	R6X39542
EyeCee One Preloaded +16.0	R7X07505
EyeCee One Preloaded +26.5	R6X06059
EyeCee One Preloaded +25.0	R6X08475
EyeCee One Preloaded +21.5	R6X12935
EyeCee One Preloaded +23.5	R6X26583
EyeCee One Preloaded +22.0	R6X34596
EyeCee One Preloaded +20.0	R6X37028
EyeCee One Preloaded +20.0	R6X37029
EyeCee One Preloaded +24.0	R6X32450
EyeCee One Preloaded +24.0	R6X32451
EyeCee One Preloaded +24.0	R6X32452
EyeCee One Preloaded +24.0	R6X32453
EyeCee One Preloaded +24.0	R6X32459
EyeCee One Preloaded +23.5	R6X37190
EyeCee One Preloaded +23.5	R6X37192
EyeCee One Preloaded +23.5	R6X37193
EyeCee One Preloaded +23.5	R6X37194
EyeCee One Preloaded +23.5	R6X37195
EyeCee One Preloaded +22.5	R7X26814

EyeCee One Preloaded +22.5	R7X26815
EyeCee One Preloaded +22.5	R7X26816
EyeCee One Preloaded +22.5	R7X26817
EyeCee One Preloaded +22.5	R7X26818
EyeCee One Preloaded +21.5	R8X00959
EyeCee One Preloaded +21.5	R8X00402
EyeCee One Preloaded +21.5	R8X00403
EyeCee One Preloaded +20.5	R8X11541
EyeCee One Preloaded +20.5	R8X11542
EyeCee One Preloaded +26.5	R8X21754
EyeCee One Preloaded +24.0	R8X26675
EyeCee One Preloaded +24.0	R8X26676
EyeCee One Preloaded +23.5	R8X34934
EyeCee One Preloaded +22.0	R9X04342
EyeCee One Preloaded +22.0	R9X36130
EyeCee One Preloaded +22.5	R0X19291
EyeCee One Preloaded +22.5	R0X19292
EyeCee One Preloaded +23.5	R9X19768
EyeCee One Preloaded +24.5	R9X26944
EyeCee One Preloaded +21.5	R0X32337
EyeCee One Preloaded +24.0	R9X01629
EyeCee One Preloaded +20.5	R0X12821
EyeCee One Preloaded +22.0	R0X19340
EyeCee One Preloaded +20.0	R0X19701
EyeCee One Preloaded +21.5	R0X19772
EyeCee One Preloaded +21.0	R0X20094
EyeCee One Preloaded +20.5	UAX08609
EyeCee One Preloaded +21.0	UAX13254
EyeCee One Preloaded +21.0	UAX13257
EyeCee One Preloaded +24.0	UAX26502
EyeCee One Preloaded +22.0	UAX28428
EyeCee One Preloaded +20.0	UAX32893
EyeCee One Preloaded +24.0	UAX02570
EyeCee One Preloaded +23.5	UAX19879
EyeCee One Preloaded +23.5	UBX39600
EyeCee One Preloaded +24.5	UCX01558
EyeCee One Preloaded +21.5	UCX13088
EyeCee One Preloaded +23.0	UCX27146
EyeCee One Preloaded +29.0	UBX10204

EyeCee One Preloaded +19.5	UCX30694
EyeCee One Preloaded +24.0	UCX32353
EyeCee One Preloaded +25.0	UEX07493
EyeCee One Preloaded +22.5	UEX10014
EyeCee One Preloaded +23.5	UEX13289
EyeCee One Preloaded +21.0	UEX15675
EyeCee One Preloaded +26.0	UCX27664
EyeCee One Preloaded +25.5	UEX17679
EyeCee One Preloaded +23.0	UEX36416
EyeCee One Preloaded +24.5	UEX28278
EyeCee One Preloaded +24.0	UEX41521