

**Wycofywane serie — zastawki hemostatyczne Access-9™ i AccessPLUS™**

Kod katalogowy	Numer serii	Kod katalogowy	Numer serii
580758001	P2073001	K09-12917	H2288582
580774001	I2219768	K09-13501AP	H2265140
580775001	I2219771	K09-YS029A	H2266470
	I2278741		H2302474
	I2226790	K09-YS040F	H2286601
IN8112/C	H2292357	K10-05457AP	H2248304
	H2287440		H2265127
	H2281404		H2311255
	H2298626	K10-05487P	H2279857
	H2305047	K10-05635	H2285675
	H2305048		H2309893
	H2309356		H2305137
H2292361	K10-05868P		H2278583
IN8152/C	H2287442	K10-05872P	H2275807
	H2281419		H2283595
	H2287443		H2309793
	H2298629		H2319811
	H2298630	K12-00264	H2285040
	H2305054	K12-01551	H2269831
	H2281405	K12-08303	H2266358
	H2305053		H2271477
	H2309362		H2295915
K04-01197	H2296388	K12-10346	H2275972
K05-00053H	H2305217	K12-10375	H2290010
K05-00109G	H2304656	K12-10724AP	H2275826
K05-00137K	H2304653		H2299830
K05-00410G	H2305216	K12-10962	H2290113
K05-00438G	H2309130	K12-11413	H2262277
K05-00502F	H2309122		H2272368
K05-00587D	H2309211	K12-YP002A	H2269321
K05-01052	H2304607		H2278279
K05-01706C	H2309133		H2283742
K05-01880K	H2306962		H2295661
K05-02130	H2304639		K12-YP030A
K05-02595	H2304614	K12-YP060A	H2272013
K05T-01863	T2313699	K12-YP061A	H2249088
K08-YP085A	H2249051	K12-YP063C	H2272382
	H2302394		H2295556
K08-YS017B	H2249052	K12-YP064A	H2288338
K08-YS055A	H2275990		H2302500
K09-09102B	H2308279	K12-YP072	H2266224
K09-10388A	H2285501	K09-12599A	H2295450

## Wycofywane serie — zastawki hemostatyczne Access-9™ i AccessPLUS™

Kod katalogowy	Numer serii	Kod katalogowy	Numer serii
K12-YP114	H2258787	MAP112/F	I2259661
	H2288909		I2279915
	H2302503		I2279894
K12-YP115	H2288352		I2176725
	H2302505		I2215955
K12-YS018A	H2284956		I2289611
K12-YS019A	H2302506		I2279895
K12-YS026A	H2276014		I2289613
	H2302507		I2289612
K12-YS047A	H2290278		MAP150/F
MAP111/F	I2259670		I2259668
	I2279884	MAP152/B	P1908669
	I2259657		P1908670
	I2289610	MAP152/F	I2215502
MAP112/B	P1908667		I2259669
			I2279933
			I2289615
			I2176727
			I2289614
			I2289616
		I2279930	